Agilent Technologies

USB/LAN/GPIB Interfaces





Contents USB/LAN/GPIB Interfaces Connectivity Guide

CONNECTIVITY GUIDELINES

Getting Started	7
Using This Guide	8
Steps to Get Started	9
Connecting Instruments to LANs	11
LAN Quick Start	12
Step 1: Select LAN Network Type 1	15
Typical Site LAN Networks 1	16
Typical Private LAN Networks 1	17
Step 2: Gather Network Information 1	18
Gather Site LAN Network Information	19
Define Private LAN Network Parameters	21
Step 3: Connect Your Instruments	23
Connect Instruments to Site LAN	24
Connect Devices to Private LAN 2	26
Step 4: Install I/O Software on Your PC 2	28
Step 5: Configure Your Instruments 2	29
Check Your Web Browser	30
Set/Change TCP/IP Parameters	30
Step 6: Configure the LAN Interface 3	33
Step 7: Communicate with Instruments 3	38
Communicating Using VISA Assistant	39
Step 8: Program Your Instruments4	10
Connecting Instruments to USB	11
USB Quick Start	12
Step 1: Install I/O Software on Your PC 4	15
Step 2: Connect Instruments to USB	16
Step 3: Assign Instrument Alias Name	18
Step 4: Check Instrument Identification	50
Step 5: Communicate with Instruments	52
Step 6: Program Your Instruments5	54



Connecting Instruments to GPIB	55
GPIB Quick Start	56
Step 1: Install I/O Software on Your PC	59
Step 2: Install GPIB Cards in Your PC	60
Step 3: Connect Instruments to GPIB Card	64
Step 4: Configure GPIB Interface Cards	68
Step 5: Communicate with Instruments	73
Step 6: Program Your Instruments	75
Programming Your Instruments	77
Programming Overview	78
Addressing Instruments	81
Addressing Instruments via LAN	81
Addressing Instruments via USB	86
Addressing Instruments via GPIB	88
	04
Example Programs	

REFERENCE INFORMATION

Installing I/O Software	
Pre-Installation Checks	
Installing Your Software	
Troubleshooting Guidelines	
Troubleshooting LAN Interfaces	
LAN Troubleshooting Overview	
Hardware Checks	111
Communication Checks	112
Web Browser Checks	114
Troubleshooting USB Interfaces	
USB Troubleshooting Overview	
USB Hardware Checks	
USB Software Checks	
Troubleshooting GPIB Interfaces	
GPIB Troubleshooting Overview	
GPIB Hardware Checks	
GPIB Software Checks	
Agilent IO Libraries Checks	



TCP/IP Network Basics	139
LAN Interface Overview	140
Typical Network Topologies	140
LAN Hardware Architecture	142
TCP/IP Protocols	144
The TCP/IP Network Model	144
The Network Interface Layer	146
The Internet Layer	147
The Transport Layer	148
The Application Layer	150
IP Addressing	153
IP Address Classes	153
Subnets and Subnet Masks	155
Local and Remote Networks	156
IP Address Configuration Methods	158
Configuration Methods Overview	158
Dynamic Host Configuration Protocol (DHCP)	159
Auto-IP/ZEROCONF	159
Duplicate IP Address Detection	159
Device Hostname Services	161
Device Hostname Services Overview	161
Dynamic DNS Naming	162
RFC NetBIOS Naming	163
Static DNS Naming	164
Configuring Your PC for LAN Operation	165
Checking PC Settings	165
Installing Network Interface Cards	166
Installing TCP/IP on Your PC	166
Setting PC IP Address	168
Setting PC Hostname	170
Guide Information	173
Guide Contents	174
Related Documentation	175
Accessing an Electronic Copy of This Guide	175
General Information for This Guide	176
Contacting Agilent	177
Glossary	179
Index	189





Getting Started

- Using This Guide
- Steps to Get Started

NOTE

This guide does not describe LAN networks that include a gateway, such as the Agilent E5810A LAN/GPIB Gateway for Windows. See the applicable gateway documentation for information on gateway systems.

Also, this guide does not provide a detailed description of LAN, USB, or GPIB interfaces or TCP/IP networks. Consult standard reference texts for this information.

If you need to contact Agilent, see <u>Contacting Agilent</u> for addresses.



Using This Guide

This Agilent Technologies USB/LAN/GPIB Interfaces Connectivity Guide shows how to connect instruments to USB, LAN, and GPIB interfaces and how to configure and troubleshoot these interfaces on PCs with Windows 98, Windows Me, Windows NT 4.0, Windows 2000, or Windows XP operating systems. A summary of the guide contents follows.

Section	Description
CONNECTIVITY GUIDELINES	
Getting Started	Shows steps to get started and how to contact Agilent
Connecting Instruments to LANs	Gives guidelines to connect instruments to LANs
Connecting Instruments to USB	Gives guidelines to connect instruments to USB
Connecting Instruments to GPIB	Gives guidelines to connect instruments to GPIB
Programming Your Instruments	Gives guidelines to program instruments via LAN/USB/GPIB
REFERENCE INFORMATION	
Installing I/O Software	Shows how to install I/O software on your PC
Troubleshooting Guidelines	Troubleshooting guidelines for LAN, USB, and GPIB Interfaces
TCP/IP Network Basics	Provides an overview of TCP/IP networks
Guide Information	Lists general information for this guide
Glossary	Defines some of the terms used in this guide
Index	Index of selected items in this guide





Steps to Get Started

This figure shows suggested steps to get started setting up LAN, USB, or GPIB interfaces and the applicable sections of this guide. Some guidelines follow this figure to select the interface you want to set up.



✓ 9

- Connect Instruments to LAN. A Local Area Network (LAN) provides instrument connectivity over distances and allows sequential sharing of instruments among multiple PCs. A LAN is the recommended method to connect instruments together in applications such as new test systems. If you select LAN connections, go to <u>Connecting</u> <u>Instruments to LANs</u>.
- Connect Instruments to USB. Universal Serial Bus (USB) is a quick and easy way to connect instruments to PCs on a benchtop. Typically, you can use USB to connect a single instrument to a PC. If you select USB connections, go to <u>Connecting Instruments to</u> <u>USB</u>.
- Connect Instruments to GPIB. The General Purpose Interface Bus (GPIB) should be used only in a system that is already using GPIB to connect instruments. GPIB is not recommended for applications such as new systems. If you select GPIB connections, go to <u>Connecting Instruments to GPIB</u>.



Connecting Instruments to LANs





LAN Quick Start

This section shows suggested steps to help you quickly get started connecting and configuring your LAN-enabled instruments for site LAN or private LAN operation.

What are Site LANs A site LAN network is defined as a local area network (LAN) in which LAN-enabled instruments and Windows 98/2000/NT 4.0/Me/XP PCs are connected to a site LAN (workgroup LAN, Intranet, or enterprise LAN) via (optional) routers, hubs, and/or switches.

A **private LAN network** is defined as a local area network (LAN) in which LAN-enabled instruments and Windows 98/2000/NT 4.0/Me/XP PCs are **NOT** connected to a site LAN. This figure shows example site LAN and private LAN networks.



Getting Started

To get started, you may want to copy the following overview of steps to use as a guide as you set up your instruments for site LAN or private LAN operation. See the associated step if you need any more information. See these sections for additional information on setting up your instruments for site LAN or private LAN operation:

- TCP/IP Network Basics for an introduction to TCP/IP networks
- Troubleshooting LAN Interfaces for LAN troubleshooting tips
- Glossary for a definition of some LAN terms
- <u>Contacting Agilent</u> if you need to contact Agilent

◀< 12 ▶</p>





🗲 🖌 13 🕨



8 Program Your Instruments

This is an optional step. As desired, you can program instruments via the LAN using applications such as IVI-COM, VXI*pnp*, or VISA.

VISA: viOpen (..."TCPIP0::instrumentHostName::INSTR"...)



Step 1: Select LAN Network Type



- What's in This Step? The first step is to determine whether you want your instruments and your PC to be connected to a site LAN or to a private LAN. This step shows some examples of site LAN and private LAN networks. Select site LAN or private LAN connections, based on your requirements. Then, go to <u>Step 2: Gather</u> <u>Network Information</u>.
 - □ For examples of typical site LAN networks, go to <u>Typical Site LAN Networks</u>
 - □ For examples of typical private LAN networks, go to <u>Typical Private LAN Networks</u>



Typical Site LAN Networks

In this guide, a **site LAN network** is defined as a local area network (LAN) in which LAN-enabled instruments and Windows 98/2000/NT 4.0/Me/XP PCs are connected to a site LAN (workgroup LAN, Intranet, or enterprise LAN) via (optional) routers, hubs, and/or switches. This figure shows some example site LAN networks.





Typical Private LAN Networks

In this guide, a **private LAN network** is defined as a local area network (LAN) in which LAN-enabled instruments and Windows 98/2000/NT 4.0/Me/XP PCs are **NOT** connected to a site LAN (workgroup LAN, Intranet, or enterprise LAN). This figure shows some example private LAN networks.

NOTE

For a private LAN you, the designer of the LAN, are the "System Administrator" for the LAN and are responsible for defining all private LAN parameters.





Step 2: Gather Network Information

2	Gather Network Information
er n bara	etwork parameters for site LAN connections or define pr meters.
	Site LAN Network Information Card
	Instrument Information (Completed by Instrument User) (Serial Number and Ethernet (MAC) Hardware Address usually on label on instrument)

What's in This Step? When you have selected the type of LAN to connect your instrument to, the next step is to gather information for your site LAN or to define parameters for your private LAN. When you have gathered the necessary information, go to <u>Step 3: Connect Your Instruments</u>.

- □ For steps to gather site LAN network information, go to <u>Gather Site LAN Network Information</u>.
- □ For steps to select private LAN network parameters, go to <u>Define Private LAN Network Parameters</u>.



Gather Site LAN Network Information

Suggested steps follow to gather information for instrument connections to a site LAN.

1 Does the Site LAN Support DHCP? Find out if the site LAN supports Dynamic Host Configuration Protocol (DHCP). In general, if the site LAN supports DHCP, you can connect your instruments to site LAN without contacting your IT department.

However, if your local policy requires it, contact your System Administrator in your Information Technology (IT) department and tell him/her you want to connect instruments to the LAN.

2 Select Automatic Configuration. If the site LAN supports DHCP, IP addresses for the instruments are automatically assigned by the network. If you do not require fixed (static) IP addresses for your instruments, you do not need to complete the *Site LAN Network Information Card*. Go to <u>Step 3: Connect Your Instruments</u>.

OR

- 2 Select Manual Configuration. If the site LAN does not support DHCP or if you want to set a fixed (static) IP address for each instrument, the System Administrator must assign the static IP addresses and you must manually configure your instruments.
- 3 Enter Instrument Information. If the LAN does not support DHCP or if you selected manual configuration, make a copy of the Site LAN Network Information Card for each instrument to be connected to the site LAN. Record the instrument serial number and hardware (MAC) address on the Instrument Information section of the card.
- 4 Give the Cards to the System Administrator. Give the partially completed card(s) to the System Administrator and ask him/her to complete the Site LAN Information part of each card and to return the card(s) to you. Then, go to <u>Step 3: Connect Your Instruments</u>.



	Site LAN Network Inf	ormation Card
Instrument Informat (Serial Number and Et Instrument Serial I Ethernet (MAC) Ha	t ion (Completed by Ins hernet (MAC) Hardware A Number: ardware Address:	a trument User) ddress usually on label on instrument)
Site LAN Informatio	n (Completed by Syste	em Administrator)
Does the Network S	Support DHCP?	Yes No
If No, provide:	IP Address (Static):	
	Subnet Mask:	<u> </u>
	Default Gateway:	
Does the Network S	Support DNS?	Yes No
If Yes, provide:	Instrument Host Nam	e:
	Instrument Domain N	ame:
	DNS Server IP Addre	SS:
Will You Allow Univ	rersal Plug&Play to be E	nabled? Yes No



Define Private LAN Network Parameters

Suggested steps follow to define parameters for a private LAN.

NOTE

For a private LAN you, the designer of the LAN, are the "System Administrator" for the LAN and are responsible for defining all private LAN parameters. However, you may want to check with the System Administrator in your Information Technology (IT) department for guidelines on designing your private LAN.

1 Does the LAN Support DHCP? In general, private LANs with a router support DHCP but typically do not support Dynamic DNS. Determine if the LAN supports DHCP.

NOTE

Most Agilent products and PCs will automatically choose an IP address via auto-IP if a DHCP server is not present.

- 2 Complete the Network Information Card. Make a separate copy of the Private LAN Network Information Card for each instrument to be connected to the private LAN. For each card, define and record the network, PC settings, and instrument settings you want for the LAN. As required, see <u>TCP/IP Network Basics</u> or the <u>Glossary</u> for a definition of terms.
- **3** Where to go Next. When you have completed the *Private LAN Network Information Cards*, go to <u>Step 3: Connect Your Instruments</u>.

4 21 **•**

	Private LAN Network Information Card
N E T W O R K	DHCP Enabled: Yes No No Dynamic DNS Enabled: Yes No Yes No UPnP Enabled OK: Yes No Subnet Mask: DNS Server IP Address:
P C	Hardware Address:
I N S T R U M E N T	Instrument Serial Number:



Step 3: Connect Your Instruments



What's in This Step? This section gives guidelines to connect your instruments and LAN interface devices such as routers, hubs, or switches to site LANs or to private LANs.

- □ To connect instruments to a site LAN, go to <u>Connect Instruments to Site LAN</u>.
- □ To connect instruments to a private LAN, go to <u>Connect Devices to Private LAN</u>.



Connect Instruments to Site LAN

This section shows typical ways you can connect LAN-enabled instruments and interface devices such as routers, hubs, and switches to a site LAN.

Two example network connections follow that you can use as templates to connect your instruments to the site LAN. Modify the connections as required for your application. When you have completed the instrument and device connections, go to <u>Step 4: Install I/O Software on Your PC</u>.

NOTE

The example connection steps assume the PC is configured for LAN operation and connected to an existing site LAN. If this is not the case and you need to configure the PC, see <u>Configuring Your PC for LAN</u> <u>Operation</u>.

Example: Direct Connection (Site LAN)

This figure shows one way to directly connect instruments to a site LAN.



- **1 Turn Power OFF.** Remove power to the instrument.
- **2 Connect Instrument to the LAN.** Use a standard CAT5 LAN cable to connect the instrument to the site LAN.
- **3 Turn Power ON.** Apply power to the instrument and verify proper instrument power-on sequence. See the instrument's *User's Guide* for information.



Example: Switch Connections (Site LAN)

This figure shows one way to connect instruments to a site LAN via an Ethernet Hub or Switch.



- **1 Turn Power OFF.** Remove power from the instruments.
- 2 Connect Instruments to Hub/Switch. If you use a hub or switch, use standard CAT5 LAN cables to connect each instrument to the hub or switch. See the hub/switch documentation for connections.
- **3 Turn Power ON.** Apply power to the hub/switch (if not already ON) and then to the instruments and verify proper power-on sequence. See the instrument's *User's Guide* for instrument information.

🗲 🖌 25 🕨

Connect Devices to Private LAN

This section shows typical ways you can connect LAN-enabled instruments and interface devices such as routers, hubs, and switches to a private LAN, assuming your PC is properly configured for LAN operation.

Two example network connections follow that you can use as templates to connect your instruments to a private LAN. Modify the connections as required for your application. When you have completed the instrument and device connections, go to <u>Step 4: Install I/O Software on Your PC</u>.

NOTE

The example connection steps assume the PC is properly configured for LAN operation. If this is not the case and you need to configure the PC, see <u>Configuring Your PC for LAN Operation</u>.

Example: Direct Connections (Private LAN)

This figure shows one way to directly connect a Windows 98/NT/2000/Me/ XP PC with a LAN-enabled instrument. Suggested connection steps follow.



- **1 Turn Power OFF.** Remove power from the instrument.
- **2 Connect Instrument to the LAN.** Use a standard CAT5 crossover cable to connect the instrument to the LAN.
- **3 Turn Power ON.** Apply power to the instrument and verify proper instrument power-on sequence. See the instrument's *User's Guide* for information on instrument power-on sequences.



Example: Switch Connections (Private LAN)

This figure shows one way to connect a Windows 98/NT/2000/Me/XP PC and LAN-enabled instruments via an Ethernet hub or switch. Suggested connection steps follow.



- **1 Turn Power OFF.** Remove power to instruments and the hub/switch.
- 2 Connect Instruments to Hub or Switch. Use standard CAT5 LAN cables to connect each instrument to the hub or switch. See the hub or switch documentation for connection information.
- **3 Turn Power ON**. Apply power to the hub/switch (if not already ON) and then to the instruments and verify proper power-on sequence. See the instrument's *User's Guide* for information on instrument power-on sequences.

4 27 **)**

Step 4: Install I/O Software on Your PC



What's in This Step? This step gives guidelines to install applicable I/O software on your PC. Go to <u>Installing I/O Software</u> for details. After the I/O software is installed, go to <u>Step 5: Configure Your Instruments.</u>



Step 5: Configure Your Instruments



What's in This Step? *This may be an optional step.* This step shows how to manually set TCP/IP parameters for LAN-enabled instruments, using the instrument's Web Pages.

- If your LAN supports DHCP and you selected automatic configuration for your instruments, skip this step and go to Step 6: Configure the LAN Interface.
- If the LAN does not support DHCP or you want to make configuration changes, you can use the steps in this section to configure your instruments. Then, go to <u>Step 6: Configure</u> <u>the LAN Interface</u>.



Check Your Web Browser

If your instrument is Web-enabled, the instrument includes a Web Server you can access using a supported Web Browser (Internet Explorer 5.0 or higher). You can then use your Web Browser and the instrument's Web Pages to view/modify network configuration parameters as required.

NOTE

If your instrument is not Web-enabled or you cannot access the instrument at its IP address, you must set the TCP/IP parameters from the instrument's front panel.

Set/Change TCP/IP Parameters

Display the Instrument Welcome Page. To display an instrument's Welcome Page, determine the instrument's IP address from the front panel display. Then, open your Web Browser. From the Web Browser address line, type 'http://<Instrument IP Address>', where <Instrument IP Address> is the IP address displayed on the front panel, and press Enter to display the instrument's Welcome Page.

For example, if the instrument's current IP address is 169.254.3.2 (as displayed on the instrument's front panel), typing http://169.254.3.2 and pressing Enter displays the instrument's Welcome Page. This figure shows a typical display for the Agilent 33220A.



Welcome Page	Web-Enabled 20) MHz Function/Arbitrary Wav	eform Generator
Browser Web Control	nformation about this We	b-Enabled Instrument	
View & Modify Configuration	Instrument:	33220A 20 MHz Function/Arbitrary Waveform Generator	
System Status	Serial Number:	US 0123456	
	Description:	Agilent 33220A (US 0123456)	
Print Display	Hostname:	a-33220a-23456.lvld.agilent.com	
Hale with	IP Address:	156.140.105.104	
this Page	VISA TCPIP Connect String:	TCPIP0::a-33220a-23456::INSTR	

2 View Current Configuration Page. To view the current configuration, click the View and Modify Configuration icon to display the Current Configuration Page. An example display for the 33220A follows.

Welcome Page	Modify Configuration
Web Control Parameter	Currently in use
View & Modify Configuration	ON
IP Address:	156.140.105.104
System Status Subnet Mask:	255.255.254.0
Print Display. Default Gateway:	156.140.104.1
DHCP Server:	156.140.104.102
DNS Server:	156.140.2.11
Hostname:	a-33220a-23456
Domain:	Md.agilent.com
Description:	Agilent 33220A (US 0123456)
Lan Keepalive Timeout:	1800



- 3 Display the Configuring the Instrument Page. From the Current Configuration Page, click the Modify Configuration box to display the Configuring the Instrument Page. You can then set/change instrument parameters as required
- 4 Set/Change TCP/IP Parameters. This figure shows a partial Configuring your 20MHz Function/Arbitrary Waveform Generator Page. For example, we will change the settings for the 33220A IP address, and default gateway from the values shown in the Configured Value column to those shown in the Edit Configuration column. To do this, we type the desired values in the Edit Configuration column and click Save to save the values. Then, we click Reboot 33220A to make the changes effective.

NOTE

Since the IP Address, Subnet Mask, and Default Gateway parameters are marked with an asterisk (*), you must click **Save** and then click **Reboot 33220A** to make the changes effective.

The **Configured Value** for the IP Address, Subnet Mask and Default Gateway are NOT necessarily the values currently set for the instrument. These values will be used during boot if DHCP is OFF or is unavailable.

Image: Sever Status Note: You must click "Save" before changes to parameters become effective. Parameters marked with an asterisk require that the instrument be rebooted before changes take effect. Image: Sever Status Undo Edits Save Reboot 33220A Factory Defaults Image: System Status Parameter Configured Value Edit Configuration Image: Print Display Parameters will be used during boot if DHCP is OFF or unavailable IP Address: * 156.140.104.232 Image: Print Display Image: Parameters will be used during boot if DHCP is OFF or unavailable IP Address: * 255.255.254.0	Welcome Page	Configuring your 20) MHz Function/Arb	itrary Waveform Generator
Browser Web Centrol Undo Edits Save Reboot 33220A Factory Defaults Image: System Status Parameter Configured Value Edit Configuration Image: System Status DHCP: * ON O OFF © ON Image: Print Display IP Address: * 156.140.104.232 156.140.104.232 Image: Print Display IP Address: * 255.255.254.0 255.255.254.0		Note: You must click "Save" before cha require that	nges to parameters become effe the instrument be rebooted befo	ctive. Parameters marked with an asterisk(*) a re changes take effect.
Vew & Modify Configuration Parameter Configured Value Edit Configuration Image: System Status DHCP: * ON O OFF Image: ON Image: Print Display IP Address: * 156.140.104.232 156.140.104.232 Image: Print Display IP Address: * 255.255.254.0 255.255.254.0	Browser Web Control	Undo Edi	ts Save Reboot 33220A	Factory Defaults
System Status DHCP: * ON OFF © ON Print Display IP Address: * 156.140.104.232 156.140.104.232 Subnet Mask: * 255.255.254.0 255.255.254.0	View & Modify Configuration	Parameter	Configured Value	Edit Configuration
Print Display These LAN parameters will be used during boot if DHCP is OFF or unavailable IP Address: * 156.140.104.232 Subnet Mask: * 255.255.254.0	System Status	DHCP: *	ON	O OFF O ON
Print Display IP Address: * 156.140.104.232 156.140.104.232 Print Display Subnet Mask: * 255.255.254.0 255.255.254.0		These LAN parameters will be used o	during boot if DHCP is OFF or unav	ailable
Subnet Mask: * 255.255.254.0 255.255.254.0	Print Display	IP Address: *	156.140.104.232	156.140.104.232
	Help with	Subnet Mask: *	255.255.254.0	255.255.254.0
Default Gateway: * 156.140.104.1 156.140.104.1		Default Gateway: *	156.140.104.1	156.140.104.1



Step 6: Configure the LAN Interface

6	Configure the LAN Interface
	Caufin to configure the TODID LAN Olicut interface. Then, add
036101	Config to configure the TCPIP LAN Client Interface. Then, add
your ins	trument identifiers to the interface software.
your ins	The TCPIP Devices dialog

What's in This Step? This step gives guidelines to use the **IO Config** utility (installed as part of the Agilent IO Libraries) to configure the LAN interface.



Steps to Configure the LAN Interface For an application that uses VISA (such as *VISA Assistant*) to communicate with instruments via the LAN, you must configure the **TCPIP LAN Client** interface by using the **IO Config** utility. Suggested steps to configure the **TCPIP LAN Client** interface follow.

> 1 Open the IO Config Utility. To configure a TCPIP LAN Client interface, click the Agilent IO Libraries Control icon (blue IO icon on the Windows taskbar) and click Run IO Config to display the IO Config main screen.

The **I/O Config** main screen allows you to select and configure ASRL (Serial), GPIB, TCP/IP, USB, and VXI interfaces.

For this example, we will configure TCP/IP interfaces.

- 🔜 Agilent IO Libraries Configuration IO Config - 🗆 × Options Help 10 Config configures and edits Aglent 10 interfaces. To configure a new interface, select the interface in "Available Interface" and cick "Configure" To edit a configured interface, select the interface in "Configured Interfaces" and cick "Edit" To automatically configure the interfaces identified with "", cick "Auto Config" Available Interface Types Configured Interfaces VISA Type Interface Description VISA Name SICL Name * Auto Config ASRL VISA LAN Client (e.g. E5810) GPIB 82350 PCL GPIB CA GPIB 82341 ISA GPIB Card OK 82357 USB to GPIB GPIB VISA LAN Client (e.g. E5810) GPIB-VXI Command Module *LAN Client (LAN Instruments) GPIB GPIB-VXI <u>H</u>elp TCPIP USB USB VISA LAN Client for USB *E8491 IEEE-1394 to VXI LAN Server (PC as Server) n/a Configure Edit... Bemove NOTE: Auto Config will configure interfaces identified with an ^{tor}
- 2 Display the TCPIP LAN Client Dialog Box. From the main screen, highlight the TCPIP LAN Client (LAN Instruments) menu item and then click <u>Configure</u> to display the LAN Client dialog box.





3 Set LAN Client Parameters. When the LAN Client dialog box appears, in almost all cases, you can accept the default settings. If you need to change any items, click the Help button for an explanation of the items. Change items as required and then click the OK button to re-display the IO Config main screen.



4 Re-Display LAN Client Dialog Box. When the IO Config main screen reappears, in the Configured Interfaces box highlight the TCP/IP device name you set in the previous step (VISA Interface Name of TCPIP0 and SICL Interface Name of Ian in this example). Then, click the <u>Edit...</u> button to re-display the LAN Client dialog box.



35

5 Display the TCPIP Devices Dialog Box. From the LAN Client dialog box, click the Edit VISA Config... button to display the TCPIP devices dialog box.

AN Client	×	
Questions? Press the Help button below. Recommended default values are shown.		
SICL Interface Name: Ian	ОК	
VISA Interface Name: TCPIP0	Cancel	
Server Timeout: 120	Help	
Client Timeout Delta: 25	Defaults	
Cog Errors	dit VISA Config	Click the Edit VISA
Default Protocol AUTO (automatically detect proto	col)	Config button
C VXI-11 (TCP/IP Instrument Protoc	col)	

6 Display the Add a TCPIP device Dialog Box. From the TCPIP devices dialog box, click the Add device button to display the Add a TCPIP device dialog box.

TCPIP devices	X	
TCPIP0 INSTR devices present:	OK	
	Add device Click the Add device button	ce
I		


7 Enter Device IP Address/Hostname. From the Add a TCPIP device dialog box, enter either the IP address or hostname of the LAN instrument to be added and then click the OK button to re-display the TCPIP devices dialog box. For example, this dialog box adds a TCP/IP device at IP address 156.140.105.104.

Add a TCPIP device	×
Enter the address of this device: Machine Name/IP Device Name TCPIP0:: 156.140.105.104 :: :INSTR TCPIP0::156.140.105.104::instr	OK Click OK
Enter IP Address	

8 Display Devices Present. The TCPIP devices dialog box displays the TCP/IP devices present on the network. For example, this dialog box shows a device with identifier TCPIP0::156.140.105.104::inst0::INSTR is present on the network. Click the OK button to close this dialog box. Then, click the OK button on each open dialog box to close the box.

TCPIP devices	×
TCPIP0 INSTR devices present: TCPIP0::156.140.105.104::inst0::INSTR	OK Cancel
	Add device Remove device

9 Repeat for Other Instruments. If you have more than one instrument to be added to the network, repeat the previous steps for each instrument to be added. Then, go to <u>Step 7: Communicate with Instruments</u>.



Step 7: Communicate with Instruments



What's in This Step? This section gives guidelines to communicate with your instruments using *VISA Assistant*. When you have verified communication with your instruments (as desired), go to <u>Step 8: Program Your Instruments</u>.

NOTE

Communication with installed LAN instruments was established in <u>Step 6: Configure the LAN Interface</u> if the instrument(s) were displayed in the **TCPIP devices** dialog box. Thus, this is an optional step you can use to verify communication with instruments.



Communicating Using VISA Assistant

VISA Assistant is an application program that uses the Agilent IO Libraries to communicate with LAN, USB, and GPIB instruments. VISA Assistant can automatically detect and assign VXI*plug&play* instrument drivers to instruments.

VISA Assistant can also be used to send and receive strings to instruments which support formatted I/O. For other instruments, VISA Assistant allows you to read and write memory areas. An example follows to show one way to communicate with a 33220A using VISA Assistant.

Example: Communicating Using VISA Assistant

This example shows one way to use *VISA Assistant* to communicate with an Agilent 33220A to return the ID string for the instrument. To open *VISA Assistant*, click the **Agilent IO Libraries Control** icon (blue icon on the Windows taskbar) and click **Run VISA Assistant** to display the *VISA Assistant* main screen. Then, use the steps in this figure to return the identification string for an Agilent 33220A.



┥ ┥ 39 🕨

Step 8: Program Your Instruments



What's in This Step? *This is an optional step.* After you have set up your site or private LAN, as desired you can program LAN-enabled instruments from your PC via the LAN. See <u>Programming Your Instruments</u> for guidelines and examples.



Connecting Instruments to USB





USB Quick Start

This section shows suggested steps to help you quickly get started connecting USB-enabled instruments to the Universal Serial Bus (USB).

NOTE

Optionally, a USB hub may be connected between the PC and USB instrument(s). However, this configuration is not described in this guide. See your USB hub documentation if you use a USB hub.

Typical USB Interface System

In this guide, a **USB interface system** is defined as a system in which USB-enabled instruments are connected via a USB cable to a USB 1.1 port in a Windows 98 (SE)/Me/2000/XP PC or to a USB 2.0 port in a Windows XP (with Service Pack 1) PC. This figure shows a typical USB interface system with a PC and a USB instrument connected via a USB cable.



Getting Started To get started, you may want to copy the figures for Steps 1 - 6 to use as a guide as you set up your USB system. See the associated step if you need more details.

Other InformationSee these sections for additional information on setting up instruments for
USB operation:

- Troubleshooting USB Interfaces for USB troubleshooting tips
- <u>Glossary</u> for a definition of some USB terms
- Contacting Agilent if you need to contact Agilent

4 4 4 4













6 Program Your Instruments

This is an optional step. As desired, you can program instruments via the USB interface using applications such as IVI-COM, VXI*pnp*, or VISA.

VISA: viOpen (...,"ArbGen",...) SICL: iopen ("ArbGen")



Step 1: Install I/O Software on Your PC



What's in This Step? As required, you should install I/O software (such as the Agilent IO Libraries), on your PC before you connect instruments to USB. See Installing I/O Software for details. Then, go to Step 2: Connect Instruments to USB.

NOTE

Version M.01.00 or later of the Agilent IO Libraries is required for USB.



Step 2: Connect Instruments to USB



What's in This Step? This step shows how to connect USB Instruments to USB ports on your PC.



Steps to Connect USB Instruments

Use the steps in this figure to directly connect USB instruments to a Windows 98(SE)/Me/2000/XP PC. If you have not yet installed I/O software on your PC, go to Installing I/O Software and install the software BEFORE you connect USB instruments to your PC. When you have made the connections for your system, go to Step 3: Set Instrument Alias Name.

NOTE

If you have not yet installed I/O software on your PC, go to <u>Step 1:</u> <u>Install I/O Software on Your PC</u> and install the software BEFORE you connect USB instruments to your PC.

Be sure to use a USB 2.0-compliant USB cable, even if you are using USB 1.1 operation.





Step 3: Assign Instrument Alias Name



What's in This Step? This step shows how to assign an Alias name for USB instruments. You can use the Alias name when addressing the instrument.



Assigning an Alias Name When you apply power to a USB instrument and then connect the USB cable between the instrument and your PC, an **Assign USB device alias** dialog box appears. From this dialog box, you can assign the instrument Alias name, as desired.

NOTE

When power is applied, a **Found New Hardware** dialog box may appear. If so, follow the on-screen instructions for installation advice. Then, return to this section.

The USB Alias name associates an alias with a specific instrument. You can use the Alias name as a VISA *rsrcName* in *viOpen()* or as a SICL *address* in *iopen()* instead of the instrument's Vendor ID, Product ID, Serial Number, and Interface Number (see <u>Step 4: Check Instrument Identification</u> for details).

Example: Assigning a USB Alias Name

This figure shows an example **Assign USG device alias** dialog box for an Agilent 33220A. To change the USB Alias name, type the new name in the **Alias name** box and then click the **OK** button. The Alias name can be used as the preferred address for the VISA *rsrcName* and the SICL address.

For this example, the Alias name has been changed to **.ArbGen.** from the default **.UsbDevice1**.. Thus, **ArbGen** can be used in viOpen() as the VISA *rsrcName* parameter, rather than **.USB0::2391::1031::MY43000029::0::INSTR**.. Or, **ArbGen** can be used in iopen() as the SICL *address* parameter, rather than the alternate address of **.usb0[2391::1031::MY43000029::0]**.

	ArbGen
	Identification: Agilent Technologies 33220A Waveform Generator
You can assign a USB Allas	Uisa Resource Name:
name from this dialog box. (If	Preferred ArbGen
you do not change the name	Alternate USB0::2391::1031::MY43000029::0::INSTR
now, you can change it later	, SICL Address String:
by using IO Config.)	Preferred ArbGen
	Alternate usb0[2391::1031::MY43000029::0]
	(OK) Cancel

49

Step 4: Check Instrument Identification



What's in This Step? *This is an optional step* that shows suggested steps follow to check instrument identification parameters. If you do not want to check instrument identification parameters, skip to <u>Step 5: Communicate with Instruments</u>.



Steps to Check1Display the IO Config Main Screen. Click the blue IO icon (on the
Windows Taskbar) and select Run IO Config to display the IO Config
main screen.Parametersmain screen.



2 Display the USB Devices Dialog Box. Highlight the USB instrument name in the Configured Interfaces box (VISA Name of USB0 and SICL Name of usb0 for the previous figure). Then, click <u>Edit...</u> to display the USB Devices dialog box.

This dialog box allows you to associate a USB Alias name with a specific instrument. The association is based on the instrument's Serial Number, Vendor ID, Product ID, and Product Identification.

Alias haine	Serial number	Vendor ID	Product ID	Identification	Add Alias
ArbGen	MY43000029	2391	1031	Agilent Technologies 3322	
					Edit Alias
					Delete Alias
					Close

3 Add/Change/Delete Aliases. You can use the USB Devices dialog box to change existing USB Alias names or to add or delete Alias names. Click the Add Alias box to display the Add Alias dialog box, click the Edit Alias box to display the Assign USB device alias dialog box, or click the Delete Alias box to delete the highlighted alias.

NOTE

- Multiple alias names may be associated with the same instrument.
- If you delete all aliases, you can restore a default alias by closing IO Config and cycling power on the instrument.



Step 5: Communicate with Instruments



What's in This Step? This is an optional step that gives guidelines to communicate with your instruments using *VISA Assistant*. When you have verified communication with your instruments (as desired), go to <u>Step 6: Program Your Instruments</u>.

NOTE

Communication with installed instruments was established in <u>Step 2: Connect Instruments to USB</u>. Thus, this is an optional step you can use to verify communication with instruments.



Example: Communicating Using VISA Assistant

VISA Assistant is an application program that uses the Agilent IO Libraries to communicate with LAN, USB, and GPIB instruments. VISA Assistant can automatically detect and assign VXI*plug&play* instrument drivers to instruments.

VISA Assistant can also be used to send and receive strings to instruments which support formatted I/O. For other instruments, VISA Assistant allows you to read and write memory areas. VISA Assistant also describes attributes that are associated with an instrument.

This example shows one way to use *VISA Assistant* to communicate with an Agilent 33220A to return the ID string for the instrument. To open *VISA Assistant*, click the **Agilent IO Libraries Control** icon (blue **IO** icon on the Windows taskbar) and click **Run VISA Assistant** to display the *VISA Assistant* main screen. Then, use the steps in this figure to return the identification string for an Agilent 33220A.





Step 6: Program Your Instruments

6 Program Your Instruments
This is an optional step. As desired, you can program instruments via the USB interface using applications such as IVI-COM, VXI <i>pnp</i> , or VISA.
VISA: viOpen (,"ArbGen",) SICL: iopen ("ArbGen")

What's in This Step? *This is an optional step.* After you have set up your USB system, as desired you can program USB instruments from your PC via the USB interface. See <u>Programming Your Instruments</u> for guidelines and examples.



Connecting Instruments to GPIB





GPIB Quick Start

This section shows suggested steps to help you quickly get started connecting GPIB instruments to the General Purpose Interface Bus (GPIB).

Typical GPIB In this guide, a **GPIB interface system** is defined as a system in which GPIB instruments are connected to a GPIB interface card in a Windows 98/2000/ NT 4.0/Me/XP PC via GPIB cables. This figure shows a typical GPIB interface system with a PC and two GPIB Instruments connected via GPIB cables.



Getting Started To get started, you may want to copy the figures for Steps 1 - 6 to use as a guide as you set up your GPIB system. See the associated step if you need more details.

Other InformationSee these sections for additional information on setting up instruments for
GPIB operation:

- Troubleshooting GPIB Interfaces for GPIB troubleshooting tips
- <u>Glossary</u> for a definition of some GPIB terms
- Contacting Agilent if you need to contact Agilent















6 Program Your Instruments This is an optional step. As desired, you can program instruments via the GPIB interface using applications such as IVI-COM, VXI*pnp*, or VISA. VISA: viopen (..."GPIB0::5::INSTR"...)



Step 1: Install I/O Software on Your PC



What's in This Step? Before you connect your instruments to GPIB, as required install I/O software, such as the Agilent IO Libraries, on your PC. See <u>Installing I/O</u> <u>Software</u> for details. Then, go to <u>Step 2: Install GPIB Cards in Your PC</u>.



Step 2: Install GPIB Cards in Your PC



What's in This Step? This step shows how to install GPIB Interface Cards (such as an Agilent 82350 PCI GPIB Interface for Windows) in your PC.

NOTE

If you have not yet installed I/O software on your PC, go to <u>Step 1: Install</u> <u>I/O Software on Your PC</u> and install the software BEFORE you install GPIB Interface Cards in your PC.



Steps to Install a GPIB Card in Your PC

Example steps follow to install an Agilent 82350 PCI GPIB Interface Card for Windows in your PC. Modify the steps as required if you install a different GPIB Interface Card in your PC. When you have installed all required GPIB Interface Cards, go to <u>Step 3: Connect Instruments to GPIB Card</u>.

CAUTION

To reduce the risk of damaging a GPIB Interface Card, only handle the card by the sheet metal frame or by its edges.

1 **Record the Card Serial Number.** Remove the card from its anti-static bag and record the Serial Number for future reference. Save the anti-static bag so you can protect the card if you need to remove the card from the PC.

For example, the 82350 Serial Number is located on the white serial number label on the card. The label contains 24 characters, with the last 8 characters (plus US) being the 82350 Serial Number. Thus, if the numbers on the serial number label are 82350-66511-4224-01-42300113, since the last 8 characters are 42300113, the 82350 Serial Number is US42300113.

2 Remove PC Cover. Remove power from the PC and from all of its peripherals. Then, remove the power cord from the PC. Unlock and remove the cover from the PC to allow access to the I/O slots. See your PC documentation for instructions.



🗲 🧹 61 🕨

3 Remove a Cover Plate. Remove one of the PC back panel cover plates. The 82350B is a 5V PCI card and will not fit in a 3.3V PCI slot or in an EISA or ISA slot. Choose a 5V PCI slot that will give adequate clearance for the GPIB connector.



Install the 82350. Insert the 82350 Interface Card edge connector into the PCI expansion slot connector of the PC. Make sure the interface is fully seated by pushing firmly on the top edge of the card with the palm of your hand. The GPIB connector should extend through the back panel opening to allow GPIB cable connection.

If you install more than one 82350, you may want to install the cards so there is at least one empty slot between every two 82350s. When inserting the 82350, be sure to hold the card by its edges. Also, be careful with the metal faceplate around the GPIB connector as the faceplate can be bent.



62

5 Replace the Cover Plate Screw. This will hold the 82350 in place. Save the blank cover plate for use if the 82350 is later removed. Replace the PC cover(s) as described in your PC documentation.





Step 3: Connect Instruments to GPIB Card



What's in This Step? This step gives guidelines to connect GPIB instruments to a GPIB Interface Card (such as an Agilent 82350) installed in your PC by using GPIB cables. When you have made the connections for your system, go to <u>Step 4: Configure GPIB Interface Cards</u>.



Steps to Connect Instruments to GPIB Cards

1

Review Connection Guidelines. The recommended method for connecting a GPIB system is linear with the system controller (PC) at one end of the system. However, a GPIB system can also be connected together in a star, linear, or a combination configuration as long as the total number of devices on the system is \leq 15 and these guidelines are followed:

- To minimize stress on connector mountings, no more than three cable connectors blocks should be stacked on top of one another. The GPIB connector screws should be fingertightened only.
- Minimize cable length as much as possible. All system devices must have tri-state drivers and must be powered on. Systems with devices not using tri-state drivers are limited to transfer rates <250 Kbytes/sec. Turning devices on or off while a system is running may cause faulty operation.
- □ For operation with data transfer rates <500 Kbytes/sec, the total length of all GPIB cables is ≤2 meters times the number of devices connected together, up to a maximum of 20 meters.</p>
- □ For operation with data transfer rates > 500 Kbytes/sec, the total length of all GPIB cables is ≤1 meter times the number of devices connected together, up to a maximum of 15 meters.
- □ The length between adjacent devices is not critical as long as the overall restriction is met. GPIB bus extenders are available that allow operation over much greater distances.
- 2 Connect GPIB Cables to the GPIB Interface Card. Connect a separate GPIB cable to each installed GPIB Interface Card using one of the following cables. Tighten the GPIB connector screws finger-tight only. (The screwdriver slots are for removal purposes only.) Two example connections follow connect a single GPIB instrument or to connect multiple GPIB instruments.



- 10833A (1 meter)
- □ 10833B (2 meters)
- □ 10833C (4 meters)
- 10833D (0.5 meter)
- □ 8120-3448 (6 meters)
- □ 8120-3449 (8 meters)
- □ Other IEEE-488 GPIB interface bus cables, as applicable

Example: Connecting a Single GPIB Instrument

This figure shows connections from a single GPIB instrument to the GPIB connector of an Agilent 82350 GPIB Interface Card installed in your PC. You may want to record the primary GPIB address of the attached instrument for future programming use. After making the connections, reconnect the PC power cord and apply power to the PC and to attached peripherals/instruments.

CAUTION

To avoid damage to the connectors, only finger-tighten the connectors.





Example: Connecting Multiple GPIB Instruments

This figure shows one way to connect three GPIB instruments to an Agilent 82350 GPIB Interface Card. You may want to record the primary GPIB address of each attached instrument for future programming use. After making the connections, reconnect the PC power cord and apply power to the PC and attached peripherals/instruments.

CAUTION

To avoid damage to the connectors, only finger-tighten the connectors.

NOTE

Although the figure shows cable connections to GPIB Instrument 1, the connection can be to any GPIB instrument in the system.



67

Step 4: Configure GPIB Interface Cards



What's in This Step? This step shows suggested actions to use the **IO Config** utility to configure a Windows 98/Me/2000/NT/XP operating system for a PC that has an Agilent 82350 GPIB Interface Card (or equivalent) installed.



Steps to Configure GPIB Interface Cards

1

Apply Power. Apply power to the PC and to the installed GPIB instruments. As Windows 98/Me/2000/NT/XP restarts, a **Found New Hardware Wizard** may start. This figure shows a typical Windows 2000 display. The display may be different for other operating systems.



2 Install Configuration Files. Use this table for the actions to take for the operating systems listed.

OS	Action
98/Me	 When the Wizard asks for the Agilent IO Libraries CD, click OK. When the Wizard asks for the hpioclas.dll file, browse to C:\windows\system\ and click OK.
2000/XP	 Click <u>Next></u> to accept the default suggestions. Click Finish to complete the installation.

Open IO Config. Click the Agilent IO Libraries Control icon (blue IO icon on the Windows taskbar) and then click Run IO Config. When the main screen appears, highlight the GPIB 82350 PCI GPIB Card menu item and then click the <u>Configure</u> button to display the 82350 PCI GPIB Card Configuration screen. (For an Agilent 82341 card, highlight GPIB 82341 ISA GPIB Card and then click <u>Configure</u>).

NOTE

For a description of IO Config and the Agilent IO Libraries, see the Agilent IO Libraries Installation and Configuration Guide. This guide may be included on the instrument CD or is available on the Web.



IO Con fig - Toci - Toci - Toei - Toar	configur onfigure dit a conf utomatica	es and edits Agilent ID interface a new interface, select the interf igured interface, select the inter ally configure the interfaces iden	s. ace in ace i ified	n 'Ava in 'Coi with "	ailable Interface nfigured Interfa '', click 'Auto Ci	Types' and click '(ces' and click 'Edit' onfig'	Configure'	
Availat	ole Interfa	ace Types		- Co	nfigured Interfa	ces		
VISA ASF ASF GPI GPI GPI GPI TCF USE USE VXI n/a	(Type 3) 3) 3) 8 8 8 8 8 8 8 8 8 8 8 9 19 19 19 19 19 19 19 19 19 19 19 19 1	Interface Description "Instance Description "Instance OPU GRID Card "Re230 CDI GRID Card Re231 ISA GRID GRID Card Re231 ISA GRID Card Card Re231 ISA GRID Card C			VISA Name	SICL Name		×Auto Config OK <u>H</u> elp
		Configure			<u>E</u> dit	Bemove	1	

4 Configure GPIB Card Parameters. When the 82350 PCI GPIB Card Configuration screen appears, set the VISA Interface Name, SICL Interface Name, Logical Unit and Bus Address values as required.

Also, verify that this is the **System Controller** for the GPIB to which it is attached (this is the typical operating mode). (See the <u>System</u> <u>Controller</u> discussion.) Then, click the **OK** button. Some guidelines to set these values follow.

Questions? Press the Help	≥ button below.
Hecommended derault value	es are snown.
Serial Number: 99307161	
SICL Interface Name: gpib0	ОК
VISA Interface Name: GPIB0 ÷	Cancel
Logical Unit: 7	Help
Bus Address: 21	Defaults
🔽 System Controller	



NOTE

After the system is configured, this screen may also display an **Edit VISA Config...** button. Clicking this button allows you to manually configure the interface as desired.

82350 GPIB Interface Card Configuration Parameters					
SICL Interface Name	Symbolic name that SICL uses to uniquely identify this GPIB interface. The default Interface Name is gpib0. The SICL Interface Name must be a unique string of alphanumeric characters, starting with a letter. Remember this value and the Logical Unit number to properly address GPIB devices in your SICL applications.				
VISA Interface Name	Symbolic name that VISA uses to uniquely identify this GPIB interface. The default VISA Interface Name is GPIB0. The 82350 Interface Name for VISA must begin with the string GPIB and have an integer appended to it, such as GPIB0, GPIB1, GPIB2, etc. Remember this value to properly address GPIB devices in your VISA applications.				
Logical Unit	Number that SICL uses to uniquely identify this 82350 interface. The Logical Unit number is an integer in the range of 0 - 10000. Remember this value and the SICL Interface Name to properly address the GPIB interface in your SICL applications.				
Bus Address	address of this GPIB interface controller on the GPIB bus. It is usually 21 if the GPIB interface is a System Controller or 20 if the GPIB interface is a non-System Controller (see System Controller, following). These addresses are chosen by convention but any address in the range 0 - 30, inclusive, may be used.				
System Controller	Determines if this interface controls which bus devices talk and which bus devices listen. If several devices exist on a bus, be sure each has a unique GPIB bus address and only one device is the System Controller (it is usually the device installed in the computer). Each GPIB interface has its own independent bus. Thus, each interface may be a System Controller as long as it is not chained together with other GPIB interfaces. However, two or more System Controllers on the same bus will cause the bus to be inoperative.				



5 Change/Accept the Configuration Values. If the configuration values displayed are acceptable to you, click the OK button. Otherwise, you can change the configuration values by clicking the arrows next to the values. If there are no arrows, you can change the configuration values by typing in the values you want.

At any time, you can click the **Defaults** button to return the dialog box to its default configuration values for the GPIB interface. When you have changed the desired values, either click the **OK** button to accept the changes or click the **Cancel** button to cancel the changes and return to the previous configuration values for the GPIB interface.

6 Repeat Steps for Other Cards. If you have installed more than one GPIB Interface Card in your PC, repeat these steps for the remaining cards. Then, go to <u>Step 5: Communicate with Instruments</u>.


Step 5: Communicate with Instruments



What's in This Step? This step gives guidelines to communicate with your instruments using VISA Assistant. When you have verified communication with your instruments (as desired), go to <u>Step 6: Program Your Instruments</u>.



Example: Communicating Using VISA Assistant

VISA Assistant is an application program that uses the Agilent IO Libraries to communicate with LAN, USB, and GPIB instruments. *VISA Assistant* can automatically detect and assign VXI*plug&play* instrument drivers to instruments.

VISA Assistant can also be used to send and receive strings to instruments which support formatted I/O. For other instruments, *VISA Assistant* allows you to read and write memory areas. *VISA Assistant* also describes attributes that are associated with an instrument.

This example shows one way to use *VISA Assistant* to communicate with an Agilent 33220A to return the ID string for the instrument. To open *VISA Assistant*, click the **Agilent IO Libraries Control** icon (blue **IO** icon on the Windows taskbar) and click **Run VISA Assistant** to display the *VISA Assistant* main screen. Then, use the steps in this figure to return the identification string for an Agilent 33220A.





Step 6: Program Your Instruments



What's in This Step? *This is an optional step*. After you have set up your GPIB system, as desired you can program GPIB instruments from your PC via the GPIB interface. See <u>Programming Your Instruments</u> for guidelines and examples.





Programming Your Instruments

- Programming Overview
- Addressing Instruments
- Example Programs



Programming Overview

This section provides an overview of programming instruments via LAN, USB, or GPIB interfaces. Four example programs are included to demonstrate generating a simple sine wave on an Agilent 33220A 20 MHz Function/Arbitrary Waveform Generator using Standard Commands for Programmable Instruments (SCPI).

The example programs are written in Microsoft Visual Basic 6.0 or Microsoft Visual C++ for Agilent VISA-COM, IVI-COM, VISA, and VXI*pnp* I/O. They are primarily intended to display SCPI commands for the 33220A and not to demonstrate the programming environment or the I/O.

NOTE

The programs in this guide are protected by copyright.

Copyright © 2003 Agilent Technologies, Inc.

You have a royalty-free right to use, modify, reproduce and distribute the programs listed in this guide (and/or any modified version) in any way you find useful, provided that you agree that Agilent has no warranty, obligations, or liability for any program content.

Agilent Technologies provides programming examples for illustration only. All example programs assume you are familiar with the programming language being demonstrated and the tools used to create and debug procedures. Agilent support engineers can help explain the functionality of Agilent software components and associated commands, but they will not modify these samples to provide added functionality or construct procedures to meet your specific needs.



Obtaining Instrument Drivers Selected combinations of program environments and I/O are included in the program examples in this guide. You will need to acquire the programming environment and the I/O independent of the *Agilent IntuiLink for the 33220A Waveform Generator* CD.

- Although not included in this guide, IDSK drivers are available and are installed as part of the free IntuiLink for Arb software included on this CD or may be obtained at http://www.agilent.com/find/intuilink.
- VISA-COM, IVI-COM, VISA, and VXIplug&play drivers for the Agilent 33220A may be obtained at http://www.agilent.com/find/33220A in the Library under the heading Software, Firmware, & Drivers.
- USB ports ONLY work properly on Windows XP when the latest Microsoft USB drivers from Windows XP Service Pack 1 have been properly installed. In addition, if using USB on Windows 2000, Agilent recommends you install the latest Service Pack for Windows 2000.

NOTE

Upgrading your XP PC to SP1 will probably not automatically update the USB drivers. See <u>www.agilent.com/find/connectivity</u> to verify if correct USB drivers are installed and to install correct drivers as required.

Installing Agilent IO Libraries

If you want to modify the example programs or write your own programs and compile them, a developer version of the Agilent E2094 I/O Libraries software must be installed. For GPIB, the Agilent E2094 I/O Libraries software is provided with Agilent GPIB I/O products. The software should have been loaded when you installed your GPIB interface card in your PC. For USB or LAN, if you do not have a copy, you can order the Agilent E2094 I/O Libraries software on CD-ROM. You must have version M (order Agilent product number E2094M) or a later version to support USB.

NOTE

You can obtain the Agilent I/O Libraries through the Agilent Developer Network. Go to **www.agilent.com/find/buyadn** and purchase an ADN Professional Membership. This gives you the right to download the latest version of the Agilent I/O Libraries. Look for the Agilent I/O Libraries link under "Downloads" on the ADN web site.



Other Visual Basic Projects

Using I/O Objects in All Visual Basic programs in this guide use the Agilent VISA-COM object. To use the I/O object in another Visual Basic project, use these steps:

- 1 Set the Reference. Set the reference to include the libraries in the Project/References menu:
 - □ "VISA COM 1.0 Type Library", corresponds to VISACOM.tlb
 - □ "Agilent VISA COM Resource Manager 1.0", corresponds to AgtRM.DLL
 - □ "VISA COM 488.2 Formatted I/O 1.0", corresponds to the BasicFormattedIO.dll
- 2 Create Formatted I/O Reference. Create the formatted I/O reference with a statement such as "Dim Fgen As VisaComLib.FormattedIO488"
- 3 Create the New Object. Use "Set Fgen = New VisaComLib.FormattedIO488" to create the actual object.



Addressing Instruments

This section gives guidelines to address instruments via a LAN, USB, or GPIB interface, including:

- Addressing Instruments via LAN
- Addressing Instruments via USB
- Addressing Instruments via GPIB

Addressing Instruments via LAN

To address instruments via the LAN, you must first configure the TCPIP LAN Client interface using the IO Config utility. A summary of the applicable LAN networking protocols follows. See <u>LAN Interface Overview</u> for a description of LAN hardware, software, and network protocols.

LAN Networking Protocols There are two LAN networking protocols provided with the Agilent IO Libraries software. You can use one or both of these protocols when configuring your systems to use VISA and SICL over LAN.

- SICL-LAN Protocol is a networking protocol developed by Agilent that is compatible with all VISA LAN products. This LAN networking protocol is the default choice in the Agilent IO Libraries configuration when configuring the LAN client. The SICL-LAN protocol on Windows 98/Me/2000/NT/XP Professional supports VISA operations over LAN to GPIB interfaces.
- VXI-11 (TCP/IP Instrument Protocol) is a networking protocol developed by the VXIbus Consortium based on the SICL-LAN Protocol that permits interoperability of LAN software from different vendors who meet the VXIbus Consortium standards.

When using either of these networking protocols, the LAN software uses the TCP/IP protocol suite to pass messages between the LAN client and the LAN server. The server accepts device I/O requests over the network from the client and then proceeds to execute those I/O requests on a local interface (such as LAN).



By default, the LAN Client supports both protocols by automatically detecting the protocol the server is using. When a VISA viOpen or SICL iopen call is performed, the LAN Client driver first tries to connect using the SICL-LAN protocol. If that fails, the driver will try to connect using the VXI-11 protocol.

If you want to control the protocol used, you can configure more than one LAN Client interface and set each interface to a different protocol. The protocol used will then depend on the interface you are connecting through.

In VISA, the protocol used is determined by the configuration settings and cannot be changed programatically. In SICL, the programmer can override the configuration settings by specifying the protocol in the *iopen* string.

Example: Using Network Protocols

Some examples using VXI-11 and SICL-LAN protocol follow. Note that SICL names are case-sensitive.

- iopen("lan[machineName]:gpib0,1") will use the configured default protocol. If AUTO is configured, SICL-LAN protocol will be attempted. If that is not supported, VXI-11 protocol will be used.
- iopen("lan; auto [machineName]:gpib0,1") will automatically select the protocol (SICL-LAN if available and VXI-11 otherwise.)
- iopen("lan;sicl-lan[machineName]:gpib0,1") will use SICL-LAN protocol.
- iopen("lan;vxi-11[machineName]:gpib0,1") will use
 VXI-11 protocol.
- The TCPIP LAN Client also supports TCP/IP socket reads and writes. To open a socket session, use iopen("lan, socketNbr[machineName]"). For example, iopen("lan, 7777 [machineName]") will open a socket connection for socket number 7777 on 'machineName'.



Configuring TCPIP LAN Client Interfaces

When you have configured VISA LAN Client interfaces, you can then use the interface name specified during configuration in a VISA viOpen call of your program. A summary of the steps to configure a TCPIP LAN Client interface follows.

NOTE

A single **TCPIP LAN Client** interface can be used by any number of **VISA LAN Client** interfaces. The main reason for having more than one LAN Client configured is that you want to use different protocols (such as SICL-LAN or VXI-11) in different **VISA LAN Clients**.

1 To configure TCPIP LAN Client interfaces using the IO Config utility, click the **Agilent IO Libraries Control** icon (blue IO icon on the Windows taskbar) and click **Run IO Config**. When the main screen appears, highlight the **TCPIP** LAN Client (LAN Instruments) menu item and then click the <u>Configure</u> button.

IO Config configure - To configure - To edit a co - To automati	ures and edits Agilent ID interfaces. e a new interface, select the interfac nfigured interface, select the interfac cally configure the interfaces identifie	e in 'Ar e in 'C ed with	vailable Interface ionfigured Interfac '*', click 'Auto Co	Types' and click 'C ces' and click 'Edit' onfig'	Configure'
- Available Inter	rface Types		Configured Interfa	ces	
VISA Type ASRL ASRL GPIB GPIB GPIB GPIB GPIB MAI TCPIP USB USB VXI N/a	Interface Description *RS-232 COM Ports VISA LAN Client (e.g. E5810) *82350 PCI GPI8 Card 82341 ISA GPI8 Card 82357 USB to GPI8 VISA LAN Client (e.g. E5810) GPI8-V4C Command Module *LAN Client (LAN Instruments) *USB Instruments VISA LAN Client (JUSB *E8491 IEEE-1394 to VXI LAN Server (PC as Server)		VISA Name	SICL Name	* <u>Auto Confi</u> OK <u>H</u> elp
	Configure		Edit	Remove	1

When the LAN Client screen appears, set the SICL Interface Name, VISA Interface Name, Logical Unit number, Server Timeout value, and Client Timeout Delta as required. Also, you can select the Default Protocol (AUTO, VXI-11, or SICL-LAN) and choose whether or not to Log Errors. Then, click the OK button.



NOTE

After the system is configured, this screen may also display an Edit VISA Config... button. Clicking this button allows you to manually configure the interface as desired.

AN Llient	×			
Questions? Press the Help button belo Recommended default values are show	ow. Wh.			
SICL Interface Name: Ian	ОК			
VISA Interface Name: TCPIP0	Cancel			
Logical Unit: 30 📩				
Server Timeout: 120	Help			
Client Timeout Delta: 25	Defaults			
V Log Errors				
Default Protocol				
 AUTO (automatically detect protocol) 				
O VXI-11 (TCP/IP Instrument Protocol)				
C SICL-LAN				

Example: Addressing LAN Instruments Using VISA/SICL

The TCPIP LAN Client interface system in this figure consists of a Windows PC with a LAN (NIC) card and three LAN instruments. Instrument1 and instrument2 are VXI-11.2 (GPIB Emulation) instruments and instrument3 is a VXI-11.3 LAN instrument.

For this system, the IO Config utility has been used to assign the LAN card a VISA name of "TCPIP0" and SICL name of "lan". For the addressing examples, instrument1 has been addressed by its machine name, instrument 2 has been addressed by its IP address, and instrument3 by its LAN name (inst0).

Since unique names have been assigned by IO Config, you can now use the VISA viOpen command or the SICL iopen command to open the I/O paths to the GPIB instruments as shown in the figure.





📢 ┥ 85 🕽

Addressing Instruments via USB

As desired, you can use an **Alias name** to address instruments via USB. The USB Alias name associates an alias with a specific instrument. You can use as a VISA *rsrcName* or SICL *address* instead of using the instrument's Vendor ID, Product ID, Serial Number, and Interface Number. See the following steps to set, add, delete or change an Alias name.

1 Display the IO Config Main Screen. Click the blue IO icon (on the Windows Taskbar) and select Run IO Config to display the IO Config main screen.



2 Display the USB Devices Dialog Box. Highlight the USB instrument name in the Configured Interfaces box (VISA Name of USB0 and SICL Name of usb0 for the previous figure). Then, click <u>Edit...</u> to display the USB Devices dialog box.



3 Add/Change/Delete Aliases. You can use the USB Devices dialog box to change existing USB Alias names or to add or delete Alias names.



Example: Assigning a USB Alias Name

This figure shows an example **USB Devices** dialog box for an Agilent 33220A. To change the USB Alias name, type the new name in the **Alias name** box and then click the **OK** button. The Alias name can be used as the preferred address for the VISA *rsrcName* and the SICL address.

For this example, the Alias name has been changed to **.ArbGen.** from the default **.UsbDevice1**.. Thus, **ArbGen** can be used in viOpen() as the VISA *rsrcName* parameter, rather than **.USB0::2391::1031::MY43000029::0::INSTR.**. Or, **ArbGen** can be used in iopen() as the SICL *address* parameter, rather than the alternate address of **.usb0[2391::1031::MY43000029::0]**..

Alias name	Serial number	Vendor ID	Product ID	Identification	Add A
ArbGen	MY43000029	2391	1031	Agilent Technologies 33220A \	
					Edit A
					Delete



Addressing Instruments via GPIB

To address instruments via a GPIB interface, you must first configure the interface using the IO Config utility.

Steps to Configure aSteps to configure a GPIB interface for an Agilent 82350 PCI card follow.GPIB InterfaceModify the steps as required if you use an 82341 ISA card.

1 To configure GPIB (82350 PCI) interfaces using IO Config, click the Agilent IO Libraries Control icon (blue icon on the Windows taskbar) and click Run IO Config. When the main screen appears, highlight the GPIB *82350 PCI GPIB Card menu item (highlight GPIB 82341 GPIB ISA Card for an 82341 card) and then click the <u>Configure</u> button.

Config confi - To configu - To edit a c - To automa	gures and edits Agilent IO interfaces. Ire a new interface, select the interfac configured interface, select the interfac tically configure the interfaces identifi	e in 'Ava ce in 'Co ed with ''	ilable Interface nfigured Interfa ", click 'Auto Co	Types' and click 'Coni ces' and click 'Edit' onfig'	figure'
Available Int	erface Types	n r-Co	nfigured Interfa	ces	
VISA Type ASRL ASRL GPIB GPIB GPIB GPIB-VXI TCPIP USB USB VXI n/a	Interface Description 'R5-232 CDM Ports VISA LAN Client (e.g. E5810) '82350 PCI GPIB Card 82341 ISA GPIB Card 82357 USB to GPIB VISA LAN Client (e.g. E5810) GPIB-VXI Command Module 'LAN Client (LAN Instruments) 'USB Instruments VISA LAN Client for USB 'E8431 IEEE-1394 to VXI LAN Server (PC as Server)		VISA Name	SICL Name	<u>* Auto Config</u> ОК <u>Н</u> еір
	Configure		<u>E</u> dit	<u>R</u> emove	



- 2 When the 82350 PCI GPIB Card Configuration screen appears, set the VISA Interface Name, SICL Interface Name, Logical Unit and Bus Address values as required. Then, click the OK button.
- 3 If you have more than one 82350 card in your system, repeat for remaining cards.

NOTE

After the system is configured, this screen may also display an **Edit VISA Config...** button. Clicking this button allows you to manually configure the interface as desired.

32350 PCI GPIB Card Configuration Questions? Press the Help	button below.
Recommended default value	ues are shown.
Serial Number: 99307161	
SICL Interface Name: gpib0	ОК
VISA Interface Name: GPIB0	Cancel
Logical Unit: 7	Help
Bus Address: 21	Defaults
✓ System Controller	



Example: Addressing GPIB Instruments Using VISA/SICL

This figure shows example VISA and SICL addressing for GPIB instruments connected to a PC via a GPIB interface.



The GPIB interface system consists of a Windows PC with two 82350 GPIB cards connected to three GPIB instruments via GPIB cables. For this system, the IO Config utility has been used to assign GPIB card #1 a VISA name of "GPIB0" and a SICL name of "gpib0".

IO Config has also been used to assign GPIB card #2 a VISA name of "GPIB1" and a SICL name of "gpib1". With these names assigned to the interfaces, the VISA/SICL addressing is as shown in the figure.

Since unique names have been assigned by IO Config, you can use the VISA viOpen command to open the I/O paths to the GPIB instruments as shown in this figure. Or, you can use the SICL iopen command to open the I/O paths shown.



Example Programs

Four example programs follow to demonstrate sine waveform generation for an Agilent 33220A 20 Mhz Function/Arbitrary Waveform Generator. The example programs are essentially repeated to show the same functionality with different programming environments and I/O. See the *Agilent IntuiLink for the 33220A Waveform Generator* CD for other program examples. The example programs are:

- Example: Simple Sine Waveform (Visual Basic and VISA-COM)
- Example: Simple Sine Waveform (Visual C++ and IVI-COM)
- Example: Simple Sine Waveform (Visual C++ and VXIpnp)
- Example: Simple Sine Waveform (Visual C++ and VISA)

Example: Simple Sine Waveform (Visual Basic and VISA-COM)

This program selects the function as "sine," and then sets the frequency, amplitude and offset of the waveform for an Agilent 33220A. The program is written in Visual Basic 6.0 and uses VISA-COM.

```
Private Sub cmdSimpleSine Click()
Dim io mgr As VisaComLib.ResourceManager
Dim Fgen As VisaComLib.FormattedIO488
Set io mgr = New AgilentRMLib.SRMCls
Set Fgen = New VisaComLib.FormattedIO488
Set Fgen.IO = io_mgr.Open(txtIO.Text)
On Error GoTo MyError
' This program sets up a waveform by selecting the waveshape
' and adjusting the frequency, amplitude, and offset.
With Fgen
  .WriteString "*RST"
                                     ' Reset the function generator
                                     ' Clear errors and status
  .IO.Clear
registers
  .WriteString "FUNCtion SINusoid" ' Select waveshape
```

┥ ┥ 91 🕨

```
' Other options are SQUare, RAMP, PULSe, NOISe, DC, and USER
  .WriteString "OUTPut:LOAD 50" ' Set the load impedance in Ohms
                                   ' (50 Ohms default)
     ' May also be INFinity, as when using oscilloscope or DMM
  .WriteString "FREQuency 2500" ' Set the frequency.
  .WriteString "VOLTage 1.2"
                                  ' Set the amplitude in Vpp.
                                  ' Also see VOLTage:UNIT
  .WriteString "VOLTage: OFFSet 0.4" ' Set the offset in Volts
     ' Voltage may also be set as VOLTage:HIGH and VOLTage:LOW for
     ' low level and high level
  .WriteString "OUTPut ON" ' Turn on the instrument output
End With
Exit Sub
MyError:
  txtError = Err.Description & vbCrLf
  Resume Next
End Sub
```

Example: Simple Sine Waveform (Visual C++ and IVI-COM)

This example program is intended for use with Microsoft Visual C++ 6.0 and the Agilent IVI-COM drivers. This program shows how to download an arbitrary waveform using binary data. The program generates a damped sine wave using 16,000 points.

NOTE

"VISA COM 1.0 Type Library" corresponds to VISACOM.tlb and "Agilent 33220 IVI COM Driver 1.0" corresponds to Agilent33220.DLL. Use the statement "#import "Agilent33220.dll" no_namespace" to import the Agilent 33220 IVI-COM Driver.

Create a variable of type IAgilent33220Ptr to hold the reference to the object. To create the object, use "Fgen.CreateInstance(__uuidof(Agilent33220))".



```
#include <stdio.h>
#include <math.h>
// Specify the default address
#define DEFAULT LOGICAL ADDRESS OLESTR("GPIB0::10::INSTR")
#import "IviDriverTypeLib.dll" no namespace
#import "IviFgenTypeLib.dll" no namespace
#import "VisaCom.tlb" no namespace
#import "Agilent33220.dll" no namespace
void main()
{
  BSTR
           TxtAddress = DEFAULT LOGICAL ADDRESS;
         Ncycles = 10; /* Define number of sinewave cycles */
  long
  double Damp factor = -5; /* Define damping factor */
  double pi = 3.14159; /* Define pi */
  long i;
  IAgilent33220Ptr Fgen;
    CoInitialize (NULL);
  Fgen.CreateInstance( uuidof(Agilent33220));
  Fgen->Initialize(TxtAddress, VARIANT FALSE, VARIANT TRUE, "");
  Fqen->Utility->Reset();// Reset the function generator
  Fgen->Status->Clear();// Clear errors and status registers
  Fgen->System->TimeoutMilliseconds = 40000; // Set timeout to 40
seconds for long download strings
  Fgen->Output->Frequency = 5000;
                                          // Output frequency is 5 kHz
  Fgen->Output->Load = 50; // Output termination is 50 Ohms
  Fqen->Output->Voltage->Amplitude = 5; // Output amplitude is 5 Vpp
    // Compute waveform
    fprintf(stderr, "Computing Waveform...\n");
  // Create an array of data
  SAFEARRAY *psa;
  // Build a SafeArray of integer values.
  SAFEARRAYBOUND bounds = \{16000, 0\};
```

93

```
psa = SafeArrayCreate (VT I2, 1, &bounds);
  // Populate the safearray
  SHORT* Waveform;
  SafeArrayAccessData(psa, (void**) &Waveform);
  for (i = 1; i \le 16000; i++)
  {
    if(i < 14000)
     {
       Waveform[i-1] = (short)(sin(2.0 * pi * (double)Ncycles *
(double)i / 14000) * 2047
         * exp(Damp factor * (double)i / 14000) + 0.5);
     }
    else
     {
       Waveform[i-1] = 0;
  }
  SafeArrayUnaccessData(psa);
  // Download data points to volatile memory
  fprintf(stderr, "Downloading Arb...\n");
   Fgen->ArbitraryWaveform->SetDataDAC(&psa);
    Fgen->System->IO->WriteString("*OPC?", VARIANT TRUE);
    Fgen->System->IO->ReadString();
  fprintf(stderr, "Download Complete\n");
    // Download data points to volatile memory
   Fgen->ArbitraryWaveform->CopyData("DAMP SIN"); // Copy arb to
    non-volatile memory
   Fqen->ArbitraryWaveform->User = "DAMP SIN"; // Select the active
    arb waveform
   Fqen->Output->Function = Agilent33220OutputFunctionUser; // Output
the selected arb waveform
    Fgen->Output->State = VARIANT TRUE;// Enable Output
  Fgen->Close();
  exit(0);
}
```

94

Example: Simple Sine Waveform (Visual C++ and VXIpnp)

This example program is intended for use with Microsoft Visual C++ 6.0 and requires VISA revision 2.0 and the Agilent Technologies 33220A VXI*pnp* instrument driver installed. This program shows how to download an arbitrary waveform using binary data. The program generates a damped sine wave using 16,000 points.

```
#include <string.h>
#include <stdlib.h>
#include <stdio.h>
#include <math.h>
#include "aq33250a.h"
#define CHECK(func)\
do {\
   ViStatus s = (func); \setminus
    if (s < 0)
    { \
  fprintf(stderr, "Error: %s returned %d\n", #func, s);\
  exit(0); \setminus
   } \
} while (0)
/* Specify the default address */
#define DEFAULT LOGICAL ADDRESS "GPIB0::10::INSTR"
void main()
ł
  ViSession vi;
  ViRsrc
              TxtAddress = DEFAULT LOGICAL ADDRESS;
  ViInt16 Waveform[16000];
  long Ncycles = 10; /* Define number of sinewave cycles */
  double Damp factor = -5; /* Define damping factor */
  double pi = 3.14159; /* Define pi */
  lonq
        i;
  ViInt16 reply;
  CHECK(aq33250a init(TxtAddress, 0, 0, &vi));
  /* This funcion initializes the instrument */
  CHECK(ag33250a reset(vi));/* Reset the function generator */
  CHECK(aq33250a dcl(vi));/* Clear errors and status registers */
    CHECK(aq33250a formBord(vi, aq33250a ARB WAV FORM BOR SWAP));
```

┥ ┥ 95 🕨

```
/* Swap data bytes (send LSB first) */
   CHECK(aq33250a freg(vi, 5000));
/* Output frequency is 5 kHz */
   CHECK(aq33250a outpLoad(vi, ag33250a OUTPUT LOAD 50));
/* Output termination is 50 Ohms */
   CHECK(aq33250a volt(vi, 5));
/* Output amplitude is 5 Vpp */
  fprintf(stderr, "Computing Waveform...\n");
  /* Calculate data points */
  for (i = 1; i \le 16000; i++)
  {
     if(i < 14000)
     ł
       Waveform[i-1] = (short)(sin(2.0 * pi * (double)Ncycles *
(double)i / 14000) * 2047
         * exp(Damp factor * (double)i / 14000) + 0.5);
     }
     else
     {
       Waveform[i-1] = 0;
     }
  /* Download data points to volatile memory from array */
  fprintf(stderr, "Downloading Arb...\n");
  CHECK(ag33250a dataDacVolatile(vi, 16000, Waveform));
    CHECK(aq33250a opc Q(vi, &reply));
  fprintf(stderr, "Download Complete\n");
   CHECK(aq33250a dataCopy(vi, "DAMP SIN"));
/* Copy arb to non-volatile memory */
   CHECK(ag33250a funcUserName(vi, "DAMP SIN"));
/* Select the active arb waveform */
   CHECK(aq33250a funcShap(vi, ag33250a OUTPUT FUNC USER));
/* Output the selected arb waveform */
CHECK(aq33250a outpSetup(vi, aq33250a OUTPUT SETUP ON));
/* Enable output */
  CHECK(ag33250a close(vi));
}
```

96

Example: Simple Sine Waveform (Visual C++ and VISA)

This example program is intended for use with Microsoft Visual C++ 6.0 and requires VISA to be installed. This program uses the arbitrary waveform function to download and output a square wave pulse with a calculated rise time and fall time. The waveform consists of 4000 points downloaded to the function generator as ASCII data.

```
#include <string.h>
#include <stdlib.h>
#include <stdio.h>
#include <visa.h>
#define CHECK(func)\
do {\
    ViStatus s = (func); \setminus
    if (s < 0)
    {\
  fprintf(stderr, "Error: %s returned %d\n", #func, s);\
  exit(0); \setminus
    } \
} while (0)
/* Specify the default address */
#define DEFAULT LOGICAL ADDRESS "GPIB0::10::INSTR"
void main()
{
  ViSession viDefaultRM, Instrument;
             TxtAddress = DEFAULT LOGICAL ADDRESS;
  ViRsrc
  ViUInt32 actual;
  char SCPIcmd[10000];
  unsigned short i;
  CHECK(viOpenDefaultRM(&viDefaultRM));
  CHECK(viOpen(viDefaultRM, TxtAddress, VI NULL, VI NULL,
&Instrument));
  /* Specify long seconds timeout for waveform download */
  CHECK(viSetAttribute(Instrument, VI ATTR TMO VALUE, 40000));
  strcpy(SCPIcmd,"*RST\n"); /* Reset the function generator */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
```



```
strcpy(SCPIcmd,"*CLS\n"); /* Clear errors and status registers */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  /* Compute waveform */
  fprintf(stderr, "Computing Waveform...\n");
  strcpy(SCPIcmd, "DATA VOLATILE");
  sprintf(SCPIcmd, "%s,%3.1f", SCPIcmd, (double)(i - 1)/5);
  for(i = 6; i <= 205; i++)
    strcat(SCPIcmd, ",1"); /* Set pulse width (200 points) */
  for(i = 206; i <= 210; i++) /* Set fall time (5 points) */
    sprintf(SCPIcmd, "%s,%3.1f", SCPIcmd, (double)(210 - i)/5);
  for(i = 211; i <= 4000; i++)</pre>
    strcat(SCPIcmd, ",0"); /* Set remaining points to zero */
  strcat(SCPIcmd, "\n");
  /* Download data points to volatile memory */
  fprintf(stderr, "Downloading Arb...\n");
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  fprintf(stderr, "Download Complete\n");
  /* Set up arbitrary waveform and output */
  strcpy(SCPIcmd, "DATA:COPY PULSE, VOLATILE\n"); /* Copy arb to
non-volatile memory */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  strcpy(SCPIcmd, "FUNCtion:USER PULSE\n"); /* Select the active
arb waveform */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
```

┥ ┥ 98 🕨

```
strcpy(SCPIcmd, "FUNCtion:SHAPe USER\n"); /* Output the selected
arb waveform */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  strcpy(SCPIcmd, "OUTPut:LOAD 50\n"); /* Output termination is 50 Ohms
*/
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  strcpy(SCPIcmd, "FREQuency 5000;VOLTage 5\n"); /* Output frequency
is 5 kHz @ 5 Vpp */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  strcpy(SCPIcmd, "OUTPut ON\n"); /* Enable output */
  CHECK(viWrite(Instrument, SCPIcmd, (ViUInt32)strlen(SCPIcmd),
&actual));
  CHECK(viClose(Instrument));
  CHECK(viClose(viDefaultRM));
}
```



Installing I/O Software

Pre-Installation Checks

Installing Your Software



Pre-Installation Checks

Before you install I/O software on your PC, verify that your system meets the minimum hardware and software requirements to install and use the Agilent IO Libraries for the USB, LAN, or GPIB interface. Adding additional RAM may improve system performance. For example, this table shows minimum hardware and software requirements to install the software contained on the *Agilent IntuiLink for the 33220A Waveform Generator* CD.

Hardware	Minimum Requirements			
PC Operation/Memory	Pentium 200 MHz operation and 64 MBytes RAM			
Free Disk Space	80 MBytes			
Ethernet Connections (LAN)	10BASE-T or 100BASE-TX compliant			
PCI Bus Slot (GPIB)	At least one PCI Bus Slot (to install Agilent 82350B)			
ISA Bus Slot (GPIB)	At least one ISA Bus Slot (to install Agilent 82341A)			
USB Port (USB)	At least one USB Port			
Software	Minimum Requirements			
Operating System	Windows 98/Me/NT 4.0 (SP5 or greater)/2000/XP (LAN/GPIB) Windows 98 (SE)/Me/2000/XP (USB)			
Web Browser (Optional)	Internet Explorer 5.0 or higher. (Internet Explorer 5.5 or higher is recommended.)			
Agilent IO Libraries Installation Sizes	 Core Installation (LAN only) Run-Time Installation (adds hardware drivers) Full Installation (without manuals) Full Installation (with manuals) Add Adobe Acrobat Reader (version 5.05) 	5 Mb 9 Mb 14 Mb 34 Mb 11 Mb (addl)		



Installing Your Software

Suggested steps follow to select and install I/O software for your application.

Determine Your Application Tasks. One way to select the I/O software for your application is to use the product life cycle approach. As shown in this figure, as the product life cycle moves from R&D to Manufacturing, tasks for the product typically change from Product Design to Product Characterization to Product Test. The associated programming levels range from no programming to minimal programming to extensive programming.

Agilent *IntuiLink* is typically used for Product Design or Product Characterization, while IVI-COM and VXI*pnp* drivers are typically used for Product Characterization or Product Test. For any software application that uses VISA (such as Agilent *IntuiLink*), a version of the Agilent IO Libraries (developer or runtime) is required.



103

2 Select Application and Agilent IO Libraries Software. The next step is to select the I/O software application(s) and Agilent IO Libraries version you want to install on your PC. This table gives guidelines for using IO software and Agilent IO Libraries version (developer or runtime). See substep 3, following, to select the Agilent IO Libraries version.

Use This Software:	When You Want to:				
Instrument Web Browser (Product Design or Product Characterization)					
Provides a "soft front panel" displayed on your Web Browser that you can use to control the instrument.	Remotely communicate with instruments from your PC, but do not need to program the instruments. (The IO Libraries are not required.)				
Agilent IntuiLink (Product Design or Product Characterization)					
Provides a "soft front panel" specific to an instrument that you can use to control instruments.	Remotely communicate with instruments from your PC, but do not need to program the instruments (The IO Libraries are required.)				
IVI-COM Drivers (Product Characterization or Pro	duct Test)				
Interchangeable Virtual Instruments (IVI)-COM Drivers provide hardware-independent programming syntax for instruments that perform the same function (DMMs, Counters, Switches, etc.).	Develop new applications using Visual Basic, Visual C++, etc. and the IVI-COM drivers. In most cases, new code must be developed for programs using IVI-COM drivers.				
VXIplug&play Drivers (Product Characterization or Product Test)					
VXI <i>plug&play</i> functions provide a set of instrument- specific functions you can use for programming instruments using Visual Basic, Visual C++, etc.	Develop or modify existing applications using Visual Basic, Visual C++, etc. You can use existing code if the same driver (or update) is supported on both the old and new instrument.				
Agilent IO Libraries (Developer) (All Tasks)					
The developer version of the Agilent IO Libraries allows you to program your instruments using Visual Basic, Visual C++, etc. and to communicate with instruments via USB, LAN, or GPIB.	Use the developer version when you want to develop or modify existing applications using Visual Basic, Visual C++, etc. and VISA. Applications written with VISA for other instruments can be ported to your instruments with little or no modification.				
Agilent IO Libraries (Runtime) (All Tasks)					
The runtime version of the Agilent IO Libraries is a subset of the Developer version. With the runtime version of the libraries, you can communicate with instruments, but you cannot program the instruments using VISA.	Use the runtime version when you want to run existing applications that use VISA to communicate with instruments, but you do not need to modify or develop the instrument communication code in the application program.				



3 Select the Agilent IO Libraries Version. As noted, there are two versions of the Agilent IO Libraries: developer and runtime. Use the developer version if you want to communicate with instruments and to program instruments using VISA. Use the runtime version if you want to communicate with instruments, but do not need to program the instruments.

This table shows the version of the Agilent IO Libraries required for LAN, USB, and GPIB interfaces and typical programming methods.

Inter- face	When to Use this Interface:	Programming Method	Need Developer IO Libraries	Need Run- Time IO Libraries	No IO Libraries Needed
LAN	A LAN is the recommended method to connect instruments together in applications such as new test systems.	Sockets			
		IVI-COM			
		VXIplug&play			
		VISA			
USB	USB is a quick and easy way to connect instruments to PCs on a benchtop. Typically, you can use USB to connect instruments to a PC for benchtop R&D testing.	IVI-COM			
		VXI <i>plug&play</i>			
		VISA*			
GPIB	GPIB should be used only in a system that is already using GPIB to connect instruments. GPIB is not recommended for new systems.	IVI-COM			
		VXIplug&play			
		VISA			

* Revision M.01.00 or later of the Agilent IO Libraries is required for USB

4 Install Your Software. This table shows the media you can use to install your selected I/O software on your PC for both runtime and developer environments. Follow the instructions on the media to install the software.

NOTE

Installing Agilent IntuiLink or the IVI-COM drivers also automatically installs a runtime version of the Agilent IO Libraries. Revision M.01.00 or later of the Agilent IO Libraries is required for programming via USB.



For This Environment:	Use This Media:	That May Include:	To Get This Media:
Runtime (cannot program instruments)	Instrument CD	- Agilent <i>IntuiLink</i> - IVI-COM Drivers - VXI <i>pnp</i> Drivers - Agilent IO Libraries (runtime) - Agilent VEE evaluation copy	Ships with the Instrument. Check the instrument CD to see which items are included on the CD.
	Web	- Agilent IO Libraries (runtime)	Download the libraries from www.agilent.com/connectivity
Developer (can program instruments)	E2094M	- Agilent IO Libraries (developer)	Order the E2094M product from Agilent. See <u>Contacting Agilent</u> for addresses.

Example: Installing Software from the 33220A CD

For example, the *Agilent IntuiLink for the 33220A Waveform Generator* CD that ships with the 33220A instrument includes these items that you can install on your PC, as desired:

- IntuiLink Waveform Editor
- ActiveX Software Components
- VXI *plug&play* Driver
- IVI-COM Driver
- Agilent IO Libraries (runtime)
- Programming Examples
- Manuals
- Agilent VEE (evaluation copy)

For this CD, installing Agilent *IntuiLink* or the IVI-COM Drivers also automatically installs a **runtime** version of the Agilent IO Libraries. (If you want to install a **developer** version of the Agilent IO Libraries, you must install the libraries from the CD included in the E2094M product.



Troubleshooting Guidelines

- Troubleshooting LAN Interfaces
- Troubleshooting USB Interfaces
- Troubleshooting GPIB Interfaces



Troubleshooting LAN Interfaces

This section shows suggested troubleshooting steps for LAN interfaces, including the following items. For information on TCP/IP networks, see <u>TCP/IP Network Basics</u>.

- LAN Troubleshooting Overview
- Hardware Checks
- Communication Checks
- Web Browser Checks

LAN Troubleshooting Overview

LAN troubleshooting guidelines follow. However, the first step is to check for normal operation, as shown in <u>Checking for Normal Operation</u>.

Checking for Normal Defore beginning LAN troubleshooting steps, you should first check to see what symptoms appear when you attempt to connect to instruments via the LAN. To see if your PC has network access to an instrument, open the Web Browser on the PC and type in the instrument's IP address (or Hostname, if known) on the web browser's address line.

If the instrument's **Welcome Page** appears, this PC has network connection to the instrument. Repeat this step for each instrument on the network. An example for the Agilent 33220A follows. If normal operation cannot be established, go to <u>Troubleshooting Flowchart</u> for a suggested sequence of troubleshooting steps.


Example: Using a Web Browser

Since the 33220A is Web-enabled, to display the 33220A Welcome Page, first determine the 33220A IP address from the instrument's front panel display. Next, open your web browser and type 'http://<33220A IP Address>', where <33220A IP Address> is the IP address displayed on the front panel. Then, press Enter to display the 33220A Welcome page.

For example, if the current IP address is 169.254.3.2, typing http:// 169.254.3.2 and pressing Enter displays the 33220A Welcome page. This figure shows a portion of an example Welcome page display for the Agilent 33220A.

Web-Enabled	20 MHz Function/Arbitrary Wav	eform Generator
Browser Web Control Information about this W	eb-Enabled Instrument	
View & Modify Configuration	33220A 20 MHz Function/Arbitrary Waveform Generator	
System Status	US 0123456	
Description:	Agilent 33220A (US 0123456)	
nt Display Hostname:	a-33220a-23456.lvld.agilent.com	
IP Address:	156.140.105.104	
Page VISA TCPIP Connect String:	TCPIP0::a-33220a-23456::INSTR	



Troubleshooting Flowchart

If normal operation cannot be established, see this flowchart for suggested troubleshooting steps. We suggest you start at <u>Hardware Checks</u>.

1 Hardware Checks	2 Communication Checks	3 Web Browser Checks
Typical Causes	Typical Causes	Typical Causes
Bad cables/connections or power not ON for PC or instruments.	Improper or Incorrect IP addresses and/or Hostnames	Web Browser settings are incorrect.
Check Operating Lights on PC/Devices	Use the ping Command	Check Proxy Settings
Check Cables/ Connections/Power	Use route add (Private LAN Only)	Set Cache and Page Refresh
		Enable Javascript
After Doing These Checks:	After Doing These Checks:	After Doing These Checks:
If the cause is not identified, see Communication Checks	- If the cause is not identified, see Web Browser Checks	 If the cause is not identified or the problem cannot be fixed, contact Agilent for
If the cause is identified as an instrument hardware problem, contact Agilent.	 If the cause is identified as an instrument problem, but the problem cannot be fixed, contact Agilent. 	support.



Hardware Checks

The first step in troubleshooting LAN networks is to check the hardware setup for all devices connected to the LAN. This table shows some possible causes and corrective actions.

After you have made the hardware checks and taken corrective actions as required, again check for normal operation as shown in <u>Checking for Normal</u> <u>Operation</u>. If normal operation has not been established, go to <u>Communication Checks</u>.

Possible Cause(s)	Corrective Action(s)
PC Suspended or Network Interface Card (NIC) not operational.	Verify that the PC is not in a Suspended power state. If the lights on the NIC are ON, the NIC is probably operational.
Ethernet Hub or Switch not ON or not operational.	Verify that the hub/switch is powered ON. If the hub/switch lights are ON, the hub/ switch is probably operational.
Instruments are not ON or are not operational.	Verify that all instruments are powered ON. If the power-on light is ON, the instrument is probably operational. Also, check the NIC lights on the instruments. If the NIC lights are ON, the hardware LAN connection to the instrument is OK.
Cables are not properly connected	Check LAN cable connections between PC, hub/switch, and instruments. Tighten connectors or replace cables as required.



Communication Checks

If there are no hardware problems, but normal communication has not been established, the next step is to make communication checks using the **ping** command. In addition, for private LANs only, you may need to use the **route add** command to establish communication. If you can establish communication using these checks, but the web browser display is not correct, go to <u>Web Browser Checks</u>.

An example to verify communication using the ping command follows.

NOTE

Using the **ping**

Command

If you cannot verify communication to connected instruments using the ping command, try using **ping** to a device at a known address. If this is successful, the problem is likely with the instrument (for site LANs) or with the hub/switch or instruments (for private LANs).

Example: Verifying Communication Using the ping Command

For example, on Windows 2000 to send a *ping* command from your PC to an instrument click **Start | Programs | Accessories | Command Prompt**. When the **Command Prompt** window opens, type *ping <IP_address>*, where *<IP_address>* is the IP address of the instrument, and then press the **Enter** key. A typical display follows where **ping** was successful.

1icrosoft Windows 2000 [Version 5.00.2195] (C) Copyright 1985-2000 Microsoft Corp.	
C://ping 156.140.104.232	
Pinging 156.140.104.232 with 32 bytes of data:	
Reply from 156.140.104.232: bytes=32 time<10ms TTL=64 Reply from 156.140.104.232: bytes=32 time<10ms TTL=64 Reply from 156.140.104.232: bytes=32 time<10ms TTL=64 Reply from 156.140.104.232: bytes=32 time<10ms TTL=64	
Ping statistics for 156.140.104.232: Packets: Sent = 4, Received = 4, Lost = 0 (0% loss), Approximate round trip times in milli-seconds: Minimum = 0ms, Maximum = 0ms, Average = 0ms	
C:\>	
4	



A typical display follows where ping was not successful.

```
Command Prompt

Microsoft Windows 2000 [Version 5.00.2195]
(C) Copyright 1985-2000 Microsoft Corp.
C:\>ping 156.140.94.200

Pinging 156.140.94.200 with 32 bytes of data:
Request timed out.
Request timed out.
Request timed out.
Ping statistics for 156.140.94.200:
Packets: Sent = 4, Received = 0, Lost = 4 (100% loss),
Approximate round trip times in milli-seconds:
Minimum = 0ms, Maximum = 0ms, Average = 0ms
C:\>
```

Using the **route add** "If you are having difficulty with a private LAN (using a hub, switch, or even a simple crossover cable), a possible cause for no communication between the PC and connected instruments is that the PC has not recognized that it is connected to two different subnets. (See <u>Local and Remote Networks</u> for details.)

One way to put the PC and the instrument on the same subnet is to use **route add** *<Instrument IP Address> <PC IP Address>* from the **Command Prompt**. If this action does not correct the problem, go to <u>Web Browser</u> <u>Checks.</u>)

NOTE

If your devices are connected to site LAN, do not use the **route add** command.

Example: Using the route add Command

For example, on Windows 2000 to send a use the **route add** command, click **Start | Programs | Accessories | Command Prompt**. When the **Command Prompt** window opens, type **route add** *<Instrument IP Address> <PC IP Address>*, where *<Instrument IP_address>* is the IP address of the instrument and *<PC IP Address>* is the IP address of the instrument, and then press the **Enter** key.



For example, if the IP address of your PC is 155.139.103.152 and the instrument address is 169.254.58.10, the two devices are probably on different subnets. To add the instrument to the subnet for the PC, use one of the following:

! route is lost when the PC is rebooted route add 169.254.58.10 155.139.103.152

or

! route persists when the PC is rebooted route -p add 169.254.58.10 155.139.103.152

Web Browser Checks

If all hardware connections are OK and you have been able to establish communication between your PC and instruments, but the web browser display is incorrect or the web browser display does not appear at all, the problem could be that the web browser settings are incorrect.

Suggesed steps follow to check your web browser settings, including the following items. If doing these steps does not correct the problem, contact Agilent for support. See <u>Contacting Agilent</u> for addresses.

- Checking Proxy Settings
- Setting Cache and Page Refresh
- Enabling Javascript

Checking Proxy Settings A proxy is a service running on a computer that takes an information request from another computer and passes the request on as if the request was its own. Typically, many companies use proxies to control access from their internal network to the external Internet.

In some situations, you may not be able to communicate with the instrument if you are using a proxy server. Typically, the web page will time out because the requested address is not accessible by the network, although the correct IP address or hostname is entered into the web browser address field.

To fix this communication problem between the web browser and the instrument web server, the browser must be informed that any requests to the instrument should **not** utilize a proxy.



Example: Setting Proxy Notification (Internet Explorer 5.0 and above)

Example steps check proxy settings for Internet Explorer versions 5.0 and above follow. The steps are for Windows 2000. Modify the steps as required for other operating systems.

- 1 Click Start | Settings | Control Panel to display the Control Panel
- 2 From Control Panel, double-click Internet Options to display the Internet Properties dialog box
- 3 From the Internet Properties dialog box, select the Connections tab
- 4 Click the LAN Settings... button under Local Area Network (LAN) settings
- 5 Under Proxy server, if "Use a proxy server" is checked, click the Advanced... button
- 6 Under Exceptions, add the IP address of the instrument or a fullyqualified domain name (such as myinstr.example.com) to the list box called "Do not use proxy server for addresses beginning with:" Use a semicolon (;') to separate multiple entries.
- 7 Click the **OK** button on the bottom of the **Internet Options** dialog box to make changes active.

Setting Cache and Page Refresh For faster performance, most web browsers cache web pages. If a page is cached, an image of the web page is stored locally. When you navigate to a page that has already been viewed, the browser will load the page from its cache rather than loading it from the network.

Usually, this process works well for static web pages. However, this process may cause problems for 'dynamically' changing web pages. This problem may occur when you are navigating using the browser's forward/back/ refresh options.

Doing this usually causes the browser to first look in its cache to see if the page exists. If the page is cached, the browser displays the page from the cache instead of going to the network to update changes. Since instrument web pages are dynamic in nature, the cache may contain content that is outdated. If you see a problem where dynamic pages are not being updated automatically, the browser should be instructed to check for newer versions of a web page on every page visit.



Example: Cache and Page Refresh (Internet Explorer 5.0 and above)

Example steps for cache and page refresh for Internet Explorer versions 5.0 and above follow. The steps are for Windows 2000. Modify the steps as required for other operating systems.

- 1 Click Start | Settings | Control Panel to display the Control Panel
- 2 From Control Panel, double-click Internet Options to display the Internet Properties dialog box.
- 3 From the Internet Properties dialog box, select the General tab
- 4 Under **Temporary Internet files**, click the **Settings...** button to display the **Settings** dialog box
- 5 Under Check for newer versions of stored pages:, click the "Every visit to the page" radio button ("Automatically" should be OK for Internet Explorer 5.5 and >)
- 6 Click the **OK** button on the **Settings** page to accept the change
- 7 Click the **OK** button on the bottom of the **Internet Options** dialog box to make changes active

NOTE

If the previous steps to not resolve the problem, you may want to clear the cache. To do this, repeat steps 1 - 4. Then, in the **Temporary Internet** *files*, click the **Delete Files**... button to display the **Delete Files** dialog box. Click **OK** to delete the files and clear the cache.



Enabling Javascript Instruments that are Web-enabled may generate web pages that depend on Javascript and Frames. For proper operation, you may need to enble Javascript.

Example: Enabling Javascript (Internet Explorer 5.0 and above)

Example steps for enabling Javascript for Internet Explorer versions 5.0 and above follow. The steps are for Windows 2000. Modify the steps as required for other operating systems.

- 1 Click Start | Settings | Control Panel to display the Control Panel
- 2 From Control Panel, double-click Internet Options to display the Internet Properties dialog box.
- 3 From the Internet Properties dialog box, select the Securities tab
- 4 Click the Custom Level... button
- 5 Under Settings, scroll down to the category called Scripting
- 6 Click the **Enable** radio button under **Active Scripting**
- 7 Click the Enable radio button under Scripting of Java applets
- 8 Click the **OK** button on bottom of page to accept security settings
- 9 When the Warning! dialog box appears, click the YES button when prompted for "Are you sure you want to change the security settings for this zone?"
- **10** Click the **OK** button on the bottom of the **Internet Options** dialog box to make changes active.



Troubleshooting USB Interfaces

This section shows suggested troubleshooting steps for Universal Serial Bus (USB) interfaces, including:

- USB Troubleshooting Overview
- USB Hardware Checks
- USB Software Checks

NOTE

The troubleshooting procedures in this section are primarily oriented toward USB instruments and the Agilent IO Libraries. Consult your instrument's User's Guide for troubleshooting details for the instrument.

USB Troubleshooting Overview

This figure shows a typical USB interface with a PC and a USB instrument. A suggested troubleshooting flowchart for USB instruments and the Agilent IO Libraries follows. We suggest you start at <u>USB Hardware Checks</u> and then go to <u>USB Software Checks</u>, as required.







USB Hardware Checks

	This section gives guidelines to make hardware troubleshooting checks for connected USB instruments. We suggest you start your troubleshooting sequence by performing the following hardware checks. If performing the hardware checks do not correct the problem, go to <u>USB Software Checks</u> .
Check Cable Connections/ Damage	Check all USB cables for good connection to the USB connector on your PC, to all USB hubs (if installed) and the USB cable connections to connected USB instruments. An improperly attached USB connector can cause the interface to malfunction.
	Check all USB cables for cuts/damage and check for bent/misaligned/ crushed connector pins. Replace cables as required. If these actions do not solve the problem, try disconnecting and reconnecting (or replacing) USB cables. If the cable checks do not solve the problem, go to <u>Check PC/Instrument Functions</u>
Check PC/ Instrument Functions	Verify that all USB hubs, and all connected USB instruments are functional and are powered ON. Verify that host computer is not in a Suspended power management state. If these actions do not solve the problem, reboot the PC. If this does not solve the problem, go to <u>Check Device Manager</u> .
Check Device Manager	You can use the Windows Device Manager to reinstall the USB instrument driver, as required. For example, with Windows 2000, go to Control Panel by selecting Start Settings Control Panel . Then, select System Hardware Device Manager . From Device Manager , select the USB Test and Measurement Device and then click Properties .
	Tab to Driver and click Reinstall Driver . This will allow the Windows Plug and Play Manager to begin searching for a driver for the instrument. Since Device Manager may have disabled the instrument, click Enable to restart the instrument.
	NOTE
	Windows may place a misbehaving USB device in an "Unknown device" category. To see if this has happened, change the View in Device Manager

category. To see if this has happened, change the View in **Device Manager** by selecting **View | Show Hidden Devices**. Next, right-click the **Universal Serial Bus Controllers Unknown device** and select **Uninstall**. Then, select **Device Manager | Action | Scan** for hardware changes.



USB Software Checks

This section provides guidelines to make software troubleshooting checks for connected USB instruments. If you have not already done so, we suggest you start your troubleshooting sequence by performing the hardware checks in <u>USB Hardware Checks</u>. If performing the hardware and then the software checks does not correct the problem, contact Agilent for instructions. See <u>Contacting Agilent</u> for addresses.

Verify Agilent IO **1 Check Agilent IO Libraries Version.** If a version of the Agilent IO Libraries Installation **Check Agilent IO Libraries Version.** If a version of the Agilent IO Libraries has been installed, a blue IO icon is normally displayed on the Windows taskbar (on the lower right-hand side of the screen).



- □ If the IO icon is displayed, click the icon and click About Agilent IO Libraries Control to display the version. The version must be M.01.00 or greater.
- If the IO icon is not displayed, a version may still be installed. To check this, click Start | Programs and look for the Agilent IO Libraries program group.
- If this group is displayed, click Agilent IO Libraries | IO Control to display the IO icon. Then, click the icon and click About Agilent IO Libraries Control to display the installed version (must be M.01.00 or greater).
- □ If neither the IO icon nor the Agilent IO Libraries program group is displayed, no Agilent IO Libraries are installed. In this case, or if the installed version is not M.01.00 or greater, you must install the newer version (see substep 2 following).



- 2 Install Agilent IO Libraries (as Required). If Version M.01.00 or greater of the Agilent IO Libraries is not installed on your PC, use this substep. Otherwise, skip to <u>Verify USB Driver Installation</u>.
 - a Remove the USB cable from the USB port on the instrument.
 - **b** Uninstall the IO Libraries from the **Agilent IO Libraries Control** (blue **IO** icon)
 - c Re-attach the USB cable to the USB port on the instrument.
 - d If this does not correct the problem, go to <u>Verify USB Driver</u> <u>Installation</u>.

 Verify USB Driver
 1
 Check for USB Driver Files. Files are listed in their default directories for Windows 2000 and Windows 98 (SE)/Me/XP. Verify that the driver files are installed. If they are installed, go to Check IO Control Operation. If not, go to substep 2 following.

NOTE

USB ports ONLY operate properly on Windows XP when the latest Microsoft USB drivers from Windows XP Service Pack 1 have been properly installed. However, upgrading your XP PC to SP1 will probably not automatically update the USB drivers. See www.agilent.com/find/ connectivity to verify if correct USB drivers are installed and to install correct drivers as required.

Windows 2000		
Program Files	C:/Program Files/Agilent/IO Libraries/drivers/agusb48832.dll	
Driver Files	C:/Winnt/system32/drivers/ausbtmc.sys	
.inf Files	C:/Winnt/inf/ausbtmc.inf	
Windows 98 (SE)/Me/XP		
Program Files	C:/Program Files/Agilent/IO Libraries/drivers/agusb48831.dll	
Driver Files	C:/Windows/system32/drivers/ausbtmc.sys	
.inf Files	C:/Windows/inf/ausbtmc.inf	



2 Uninstall/Reinstall the Agilent IO Libraries. If the driver files are not found, uninstall the Agilent IO Libraries by inserting the *Instrument CD* and following the instructions to remove the libraries. Then, install Version M.01.00 or greater of the Agilent IO Libraries.

Check IO Control When the Agilent IO Libraries were installed, an **IO Control** was created. Operation When the **IO Control** is active, it is displayed as a blue **IO** icon on the Windows taskbar. If the **IO Control** is deactivated, SICL/VISA applications that are running with the instrument will be unable to open sessions.

By default, the **IO Control** is always active after the Agilent IO Libraries are installed and the blue **IO** icon is displayed. However, the **IO Control** may be active even though the blue **IO** icon is not displayed. There are two ways that the blue **IO** icon can be hidden:

- Clicking Hide Agilent IO Libraries Control. Clicking the blue
 IO icon and then clicking Hide Agilent IO Libraries Control
 hides the blue IO icon, but does not deactivate the IO Control.
- Clicking Exit. Clicking the blue IO icon and then clicking Exit causes a dialog box to appear that asks you if you want to terminate the Agilent IO Libraries Control. Clicking Yes hides the blue IO icon and deactivates the IO Control.

If the blue IO icon is not displayed, either the blue IO icon display has been turned off and/or the IO Control is not active. In this case, click **Start | Programs | Agilent IO Libraries** and then click **IO Control** to re-start the IO Control and to display the blue IO icon.



Troubleshooting GPIB Interfaces

This section shows suggested troubleshooting steps for an Agilent 82350 GPIB interface, including:

- GPIB Troubleshooting Overview
- GPIB Hardware Checks
- GPIB Software Checks
- Agilent IO Libraries Checks

NOTE

The troubleshooting procedures in this section are primarily oriented toward an Agilent 82350 GPIB card and the Agilent IO Libraries. Consult your instrument's User's Guide for troubleshooting details for the instrument.

GPIB Troubleshooting Overview

This figure shows a typical GPIB system with a PC and two GPIB instruments. A suggested troubleshooting flowchart for GPIB cards, installed instruments, and the Agilent IO Libraries follows. We suggest you start at Step 1 and then go to Step 2 and then to Step 3, as required.









GPIB Hardware Checks

This section gives guidelines to make hardware troubleshooting checks for the GPIB card (such as an Agilent 82350), including:

- Check Cables/Connections/Power
- Check Device Manager
- If Sound Card Does Not Work
- If Data Transfers to Devices Fail
- If IO Config Finds 82350 Card with Serial Number ffffffff

Check Cables/ We suggest you start your troubleshooting sequence by performing the following hardware checks. If the hardware checks do not solve the problem, see <u>GPIB Software Checks</u>.

NOTE

There are no user-servicable parts for the 82350. If you suspect a hardware failure for the 82350, contact Agilent for instructions to return the unit. See <u>Contacting Agilent</u> for telephone numbers/web site address.

- 1 Check GPIB Cable Connections. Check all GPIB cables for good connection to the GPIB connector on the 82350 installed in your PC and the GPIB cable connections between all connected GPIB instruments. An improperly attached GPIB connector can cause the bus to malfunction.
- 2 Check GPIB Cables for Damage. Check all GPIB cables for cuts/ damage and check for bent/misaligned/crushed connector pins. Replace cables as required.
- **3 Disconnect/Reconnect GPIB Cables.** If Steps 1 and 2 do not solve the problem, try disconnecting and reconnecting (or replacing) GPIB cables.
- 4 Check PC/Instrument Power-on. Verify that the PC and all connected GPIB instruments are functional and are powered ON. Verify that host computer is not in a Suspended power management state.
- 5 Reboot Your PC. If doing Steps 1, 2, 3 or 4 does not solve the problem, reboot the PC. If this does not solve the problem, go to <u>Check Device Manager</u>.

▲ 126 ►

Check DeviceYou can use the Windows Device Manager to reinstall the 82350 or
equivalent, as required. For example, with Windows 2000, go to Control
Panel by selecting Start | Settings | Control Panel.

Then, select **System | Hardware | Device Manager**. From Device Manager, select **82350** and then **Properties**. Tab to **Driver** and click **Reinstall Driver**. This will allow the Windows Plug and Play Manager to begin searching for a driver for the 82350. Since Device Manager may have disabled the 82350 device, click **Enable** to restart the 82350. If this does not resolve the problem, go to <u>GPIB Software Checks</u>.

- If your sound card stops working after configuring interfaces, disable the Not Work If your sound card stops working after configuring interfaces, disable the auto-detection process in IO Config. The auto-detection process in IO Config can cause sound and other cards to stop responding. To disable the auto-detection process, edit *c:\Program Files\agilent\IO Libraries\ intfcfg\intfcfg.txt* and follow the instructions to set all GPIB cards to be manually detected by IO Config.
- If Data Transfers to If devices attached to your ISA (82341) GPIB card are detected, but data transfers to those devices fail, this is usually caused either by an IRQ conflict with another card or by your system having no IRQs available for this card. If no IRQs are available, the IRQ value shown in IO Config will be -1.

For Windows NT/2000/XP: Try changing the IRQ setting to other values. You must reboot after changing the value to test it. If none work, you will have to free up other IRQs by disabling devices such as serial ports or sound cards.

For Windows 98/Me: You will have to free up other IRQs by disabling devices such as serial ports or sound cards. This is usually done in the BIOS setup for your computer.

If IO Config FindsIf IO Config reports finding an 82350 card with serial number ffffffff, this is82350 Card withtypically caused by getting into a situation where PCI cards are not properlySerial Number fffffffconfiguring. Try the following:

- 1 Check BIOS Setting. For Windows NT, make sure the BIOS Operating System setting is set to Not PnP OS or to Running Windows NT or to Other. It should not be set to Running Windows 95, Running Windows 98, Running Windows 2000, or PnP O/S
- 2 Upgrade your system BIOS to the latest version. Note that even brand new computers tend to have new BIOSs available.



3 If your computer locks up or freezes after installing. Typically, this is mostly interrupt conflicts with other drivers in the system. PCI allows sharing of IRQs, but this also means that the drivers for cards with which we are sharing an IRQ must be handling interrupt chaining properly.

Both NT and Windows 9x allow multiple drivers to install interrupt service routines (ISRs) for a single IRQ. The OS will keep a list of all the ISRs installed for each IRQ. When the IRQ is asserted, the OS will call each ISR in turn until one of them returns TRUE (meaning that it handled the interrupt). The ISR's responsibility is to correctly return TRUE if its device was interrupting, or FALSE if not.

- 4 **Perform Driver Work-Arounds.** Here are some ways to work around drivers that are not behaving properly:
 - Upgrade the drivers for devices sharing an IRQ with Agilent, including but not limited to, your video drivers, your LAN drivers, Agilent IDE and/or SCSI drivers, and your sound drivers.
 - Try to force the Agilent driver to be installed earlier in the ISR chain.
 - For Windows NT, Click Control Panel | Devices. Select the HP350i32. Select Startup... and set to System instead of the default Automatic. (Running IO Config will set this back to Automatic.)
 - □ For Windows 9x, change your BIOS **Installed OS** setting. Try all the different values and use the one that works.
 - □ For Windows 9x, make sure all devices sharing an IRQ with our card have drivers installed. Right-click **My Computer** and select **Properties**, then **Device Manager**. Highlight **Computer** and click **Properties**. Find the 82350 card and check that all other cards on the same IRQ have a valid driver, not the big yellow question-mark.
- 5 Re-Configure Your PC. Configure your PC so as to not share IRQ lines. Many PCI cards have bugs when sharing IRQ lines. You may or may not be able to do this on all PCs. Many PCs can be configured using the setup option when the PC is first booting.



- 6 **Upgrade your system BIOS.** Note that even new computers may have a newer BIOS available.
 - Make sure the BIOS Installed O/S setting is set correctly. This determines what software will configure all the VXIplug&play cards in your system. Either the BIOS or the O/S can perform the task of querying all the cards to determine their resource needs, picking a valid configuration for all these cards, and telling the cards what their actual resource settings are.
 - If the BIOS Installed O/S is set to Windows 95, Windows 98, Windows 2000 or Running a PnP O/S, the BIOS will not perform this task and will leave it to the O/S to do this. If the Installed O/S is set to Windows NT, NOT PnP O/S, or Other, the BIOS will perform this task.
 - For NT version 4.0 or earlier, the BIOS must perform this task as the O/S does not know how to do this. For Windows 9x, either the BIOS or the O/S can do this task, so try both. If you are running Windows NT, set to Running Windows NT, Not VXIplug&play OS or Other.



GPIB Software Checks

This section provides guidelines for 82350 software checks, including:

- Check for 82350 Driver Files
- Disable IO Config Auto-Detection
- Check BIOS/Interrupts Settings
- Set 82350 Read/Write Performance Mode

Check for 82350 After installing the Agilent IO Libraries, check for installed 82350 driver files. Driver Files

1 Check for 82350 Driver Files. Files are listed in their default directories.

Windows 2000		
Program Files	C:/Program Files/Agilent/IO Libraries/drivers/ag350i32.dll	
Driver Files	C:/Winnt/system32/drivers/agt82350.sys	
.inf Files	C:/Winnt/inf/agt350.inf	
Windows 98 (SE)/Me/XP		
Program Files	C:/Program Files/Agilent/IO Libraries/drivers/ag350i31.dll	
Driver Files	C:/Windows/system32/drivers/agt82350.sys	
.inf Files	C:/Windows/inf/agt350.inf	

2 Uninstall/Reinstall the Agilent IO Libraries. If the driver files are not found, uninstall the Agilent IO Libraries by inserting the *Instrument CD* and following the instructions to remove the libraries. Then, follow the instructions to re-install the libraries.

Disable IO Config Auto-Detection If your PC sound card stops working after configuring interfaces, disable the auto-detection process in IO Config, as the auto-detection process in IO Config can cause sound and other cards to stop responding. To disable the auto-detection process, edit *C:\Program Files\agilent\IO Libraries\ intfcfg\intfcfg.txt* and follow the instructions (comments in the *intfcfg.txt* file) to set all GPIB cards to be manually detected by IO Config.



Check BIOS/ Interrupts Settings If IO Config reports finding an 82350 card with Serial Number ffffffff, this is typically caused by PCI cards not properly configuring. Try the following steps. If these steps do not work, remove and re-install the 82350 and then reconfigure the card.

- 1 Check BIOS Setting. For Windows NT, make sure the BIOS Operating System setting is set to Not PnP OS or to Running Windows NT or to Other. It should not be set to Running Windows 95, Running Windows 98, Running Windows 2000, or PnP O/S
- 2 Upgrade your system BIOS to the latest version. New computers oftentimes have newer BIOSs available.
- If your computer locks up or freezes after installing. This is typically caused by interrupt conflicts with other drivers in the system. PCI allows sharing of IRQs, but this also means the drivers for cards which share an IRQ must be handling interrupt chaining properly. Both NT and Windows 9x allow multiple drivers to install interrupt service routines (ISRs) for a single IRQ. The OS keeps a list of all the ISRs installed for each IRQ.

When the IRQ is asserted, the OS calls each ISR in turn until one of them returns TRUE (meaning that it handled the interrupt). The ISR's responsibility is to correctly return TRUE if its device was interrupting or FALSE if not. Drivers that return TRUE, even though they did not service the interrupt, will cause problems. Try the following steps to perform driver workarounds:

- Upgrade the drivers for devices sharing an IRQ with Agilent, including, but not limited to, your video drivers, your LAN drivers, Agilent IDE and/or SCSI drivers, and your sound drivers.
- For Windows NT, Click Control Panel | Devices. Select the HP350i32. Select Startup... and set to System instead of the default Automatic. (Running IO Config will set this back to Automatic.)
- ❑ For Windows 9x, change your BIOS Installed OS setting. Try all the different values and use the one that works. For Windows 9x, make sure all devices sharing an IRQ with your card have drivers installed.



Right-click **My Computer** and select **Properties**, then **Device Manager**. Highlight **Computer** and click **Properties**. Find the 82350 card and check that all other cards on the same IRQ have a valid driver, not the big yellow question-mark.

- 4 Re-Configure Your PC. Configure your PC so as to not share IRQ lines. Many PCI cards have bugs when sharing IRQ lines. You may or may not be able to do this on all PCs. Many PCs can be configured using the setup option when the PC is first booting.
- 5 **Upgrade your system BIOS.** New computers may have a newer BIOS available. When installing a new system BIOS:
 - Make sure the BIOS Installed O/S setting is set correctly. This determines what software will configure all the VXIplug&play cards in your system. Either the BIOS or the O/S can perform the task of querying all the cards to determine their resource needs, picking a valid configuration for all these cards, and telling the cards what their actual resource settings are.
 - □ If the BIOS Installed O/S is set to Windows 95, Windows 98, Windows 2000 or Running a PnP O/S, the BIOS will not perform this task and will leave it to the O/S to do this. If the Installed O/S is set to Windows NT, NOT PnP O/S, or Other, the BIOS will perform this task and the system may not work properly.
 - For NT version 4.0 or earlier, the BIOS must perform this task as the O/S does not know how to do this. For Windows 9x, either the BIOS or the O/S can do this task, so try both. If you are running Windows NT, set to Running Windows NT, Not VXIplug&play OS or Other.

The 82350 card read and write calls use one of two modes:

Polling. Bytes are transferred to/from the card, one at a time. Polling mode is advantageous for transferring a small number of bytes because the setup overhead is very low, but it does require CPU involvement for each byte transferred.

Set 82350 Read/ Write Performance Mode



Interrupt. An entire buffer is transferred to/from the card without CPU involvement. Interrupt mode is advantageous for transferring large buffers because the higher per byte transfer rate more than compensates for the relatively long interrupt setup overhead.

The default behavior of the 82350 driver is to use **Polling** mode for transfers of 256 bytes or less and to use **Interrupt** mode for larger transfers. You can modify this default behavior by doing the following:

SICL: The SICL ihint(id, hint) function can be called to modify the read/ write behavior for on a SICL session. The hint values allowed are:

- I_HINT_DONTCARE (default value) Use Interrupt mode for transfer requests larger than 256 bytes, otherwise, use Polling mode.
- I_HINT_USEPOLL Use the Polling mode.
- I HINT IO Use the Interrupt mode.

VISA: The VISA visetAttribute(vi, VI_ATTR_DMA_ALLOW_EN, attrValue) function can be called to modify the read/write behavior for a VISA session. The VI_ATTR_DMA_ALLOW_EN values allowed are:

- VI_TRUE (default value) Use Interrupt mode for transfer requests larger than 256 bytes, otherwise, use Polling mode.
- VI FALSE Use the Polling mode.

Some additional factors to consider are:

- The settings discussed above are per session. This means you can open multiple sessions to a device and set different transfer modes for different sessions. The actual mode used will then depend on which session you are using for the read/write calls.
- In both SICL (with hint = I_HINT_DONTCARE) and VISA (with VI_ATTR_DMA_ALLOW_EN = VI_TRUE), the size of the read request (as specified by bufsize in a SICL iread() or count in a VISA viRead() function call) will determine the mode used even if the number of bytes actually read is less.
- The default formatted IO read buffer size is 4096 so when using this default size, formatted reads in SICL (with hint = I_HINT_DONTCARE) and VISA (with VI_ATTR_DMA_ALLOW_EN = VI_TRUE) will use Interrupt mode even when a small number of bytes are expected.



- The default formatted IO write buffer size is 128 so when using this default size, formatted writes in SICL (with hint = I_HINT_DONTCARE) and VISA (with VI_ATTR_DMA_ALLOW_EN = VI_TRUE) will used Polling mode even when a large number of bytes are being sent.
- In SICL, Polling mode will always be used for the iread(), ifread() and iscanf() regardless of the above settings, when a termchr is set (itermchr() is not set to -1).
- In VISA, Polling mode will always be used for viRead(), viBufRead() and viScanf() regardless of the above settings, when VI_ATTR_TERM_CHAR_EN = VI_TRUE.

The **crossover point** at which the Interrupt mode becomes faster then the Polling mode depends on the CPU speed, with a faster CPU having a higher crossover point.



Agilent IO Libraries Checks

This section gives guidelines to make troubleshooting checks for the Agilent IO Libraries, including:

- Check IO Libraries Installation
- Check IO Control Operation
- Install IO Libraries (if 82350 Was Installed First)

Check IO Libraries Installation Start your Agilent IO Libraries troubleshooting sequence by verifying IO Libraries installation. If the IO Libraries are installed, go to <u>Check IO Control</u> <u>Operation</u>.

> 1 Check Agilent IO Libraries Version. If a version of the Agilent IO Libraries has been installed, a blue IO icon is normally displayed on the Windows taskbar (on the lower right-hand side of the screen).



- If the IO icon is displayed, click the icon and click <u>About</u> Agilent IO Libraries Control to display the version. The version must be M.01.00 or greater.
- □ If the IO icon is not displayed, a version may still be installed. To check this, click Start | Programs and look for the *Agilent IO Libraries* program group.
- If this group is displayed, click Agilent IO Libraries | IO Control to display the IO icon. Then, click the icon and click <u>About</u> Agilent IO Libraries Control to display the installed version (must be M.01.00 or greater).
- □ If neither the IO icon nor the Agilent IO Libraries program group is displayed, no Agilent IO Libraries are installed. In this case, or if the installed version is not M.01.00 or greater, you must install the newer version (see Step 2, following).
- 2 Install Agilent IO Libraries (as Required). If Version M.01.00 or greater of the Agilent IO Libraries is not installed on your PC, install the libraries. Otherwise, go to <u>Check IO Control Operation</u>.



Check IO Control Operation When the Agilent IO Libraries were installed, an IO Control was created. When the IO Control is active, it is displayed as a blue IO icon on the Windows taskbar. If the IO Control is deactivated, SICL/VISA applications that are running with the 82350 will be unable to open sessions.

By default, the IO Control is always active after the Agilent IO Libraries are installed and the blue IO icon is displayed. However, the IO Control may be active even though the blue IO icon is not displayed. There are two ways that the blue IO icon can be hidden:

- Clicking Hide Agilent IO Libraries Control. Clicking the blue IO icon and then clicking Hide Agilent IO Libraries Control hides the IO icon, but does not deactivate the IO Control.
- □ Clicking Exit. Clicking the blue IO icon and then clicking Exit causes a dialog box to appear that asks you if you want to terminate the Agilent IO Libraries Control. Clicking Yes hides the blue IO icon and deactivates the IO Control.

If the blue IO icon is not displayed, either the blue IO icon display has been turned off and/or the IO Control (and associated iprocsvr.exe) is not active. In this case, click Start | Programs | Agilent IO Libraries and then click IO Control to re-start the IO Control and to display the blue IO icon.

Install IO Libraries (if 82350 was Installed First) If you installed the 82350 before installing the Agilent IO Libraries software, use the applicable step for your operating system to install the Agilent IO Libraries. After the drivers for all new cards are installed, you can run **Setup** from the *Instrument* CD to install and configure the Agilent IO Libraries.

NOTE

On older operating systems, you may be asked insert the 'HP I/O Libraries' CD. You should use the CD that contains the Agilent IO Libraries in this case. The card will be identified as a Hewlett-Packard card. This is necessary for backward compatibility.



- □ Windows 2000: The Found New Hardware Wizard will identify this card as a PCI Simple Communications Controller. The driver is typically located in the \Windows 2000 directory on the Instrument CD.
- ❑ Windows NT 4.0: Since Windows NT is not a VXI*plug&play* OS, all necessary driver installation and configuration is performed when you install the Agilent IO Libraries.
- □ Windows 98: The Add New Hardware Wizard will identify this card as a PCI Communication Device. The files needed for this device are typically located at the root of the *Instrument* CD.





TCP/IP Network Basics

- LAN Interface Overview
- <u>TCP/IP Protocols</u>
- IP Addressing
- IP Address Configuration Methods
- <u>Device Hostname Services</u>
- <u>Configuring Your PC for LAN Operation</u>

NOTE

The information in this section is a summary of TCP/IP networks and LANs and is not intended to be a complete discussion of the subject. Consult standard reference texts for further details on TCP/IP and LANs.



LAN Interface Overview

This section provides an overview of Local Area Networks (LANs) that use Transmission Control Protocol/Internet Protocol (TCP/IP), including:

- Typical Network Topologies
- LAN Hardware Architecture

Typical Network Topologies

LANs using TCP/IP can be divided into two categories: Site LANs and Private LANs.

Site LAN Topology A Site LAN can be a workgroup LAN, Intranet, or enterprise (corporate) LAN. Typically, a site LAN includes several PCs, routers, and servers. The LAN network is usually administered and controlled by a **System Administrator** (or **network administrator**) who installs, monitors, and troubleshoots the network. This figure shows a portion of a typical site LAN.





Within the site LAN, every device (PC, router, server, etc.) is called a **host**. A host is any device on a network that has a TCP/IP address. TCP/IP addresses are called **IP addresses**, and each device on the network must have a unique IP address. A typical IP address for a device is 156.140.105.50. In general, communication within the Site LAN and to Private LANs is behind a **firewall**.

Communication among devices on the Site LAN is controlled by **routers**. A router is a host that interfaces with other networks and can move data (called packets) from one network to another or can move data within the network, depending on the IP addresses of the hosts that are communicating with each other.

The servers on the Site LAN can supply addressing and naming services, such as Dynamic Host Configuration Protocol (DHCP) to automatically assign IP addresses to hosts on the network or Domain Name Service (DNS) or Windows Internet Naming Service (WINS) that allow hostnames to be associated with their IP addresses.

Private LAN A Private LAN (also called a local or isolated network) can operate as a local network (not connected to any other network) or can be connected to a site LAN or to another local network. This figure shows an example private LAN with three hosts (PC and two instruments). Communication between the PC and the instruments is via a **switch** or **hub**. You can also connect a single instrument to a PC by using a **crossover cable**.



M 🖌 141

LAN Hardware Architecture

A LAN is a way to extend the control of instrumentation beyond the limits of typical instrument interfaces. You can communicate with instruments using a web browser, the Telnet utility, or sockets. However, to program (send SCPI commands) over the LAN, you must first configure the LAN interface using the Agilent **IO Config** utility provided with the Agilent IO Libraries.

Client-Server Model The LAN software uses the **client-server model** of computing. Client-server computing refers to a model where an application (the **client**) does not perform all necessary tasks of the application itself. Instead, the client makes requests of another device (the **server**) for certain services

As shown in this figure, a LAN client, such as a PC, makes requests over the network to a LAN server, such as a PC or a LAN-enabled instrument. For example, in this figure the Server PC and the instruments can act as LAN servers. Once the LAN server has completed the requested operation on the instrument or device, the LAN server sends a reply to the LAN client. This reply contains requested data and status information that indicates whether or not the operation was successful.





Packet Switching A TCP/IP network is a **packet-switched** network. This figure shows an example packet-switched network. In this type of network, the computer that is sending the data (source host) breaks the data into smaller segments, called **packets**.

Each packet is individually addressed and is sent to the destination (destination host.) The destination host then reassembles the packets into the original message. Each packet can be sent from source to destination in any of several routes. The **routers** in the network control the paths of the packets.





TCP/IP Protocols

This section summarizes protocols for LAN networks that use TCP/IP for communication between hosts (devices such as computers, printers, and instruments), including:

- The TCP/IP Network Model
- The Network Interface Layer
- The Internet Layer
- The Transport Layer
- The Application Layer

The TCP/IP Network Model

The TCP/IP network model is based on **protocols** and an associated set of layers that control the actions of the network. A **protocol** is a rule or set of rules and standards for communication and data transfer between **hosts** on a network. A **host** is any device on the network (such as a computer, server, printer, etc.) that has a TCP/IP address. To exchange data on the network, the source and destination hosts must agree on the protocol and each host must have a unique (TCP/IP) address. When protocols are grouped together, they form a **protocol suite** and work together as a **protocol stack**.

TCP/IP Layers

The TCP/IP network model is based on the Department of Defense (DoD) model. This figure summarizes the functions of the layers in the model.




TCP/IP Protocol
OverviewThis table summarizes the components of each layer and shows the
applicable Request for Comments (RFC) for each component.

Layer	Component Name	Description	Specification
Internet	IP - Internet Protocol	The lowest level protocol by which data are sent from one computer to another on the Internet.	RFC 791* RFC 2460**
	ICMP - Internet Control Message Protocol	A message control and error-reporting protocol between IP devices on a network.	RFC 792
	ARP - Address Resolution Protocol	A protocol for translating IP addresses into Ethernet addresses.	RFC 826
Transport	TCP - Transmission Control Protocol	A connection-oriented protocol running over IP that reliably delivers data.	RFC 793
	UDP - User Datagram Protocol	A datagram-oriented protocol running over IP	RFC 768
Application	DNS - Domain Name System	A protocol for translating a Hostname into an IP address.	RFC 1034 and RFC 1035
	DHCP - Dynamic Host Configuration Protocol	A protocol for obtaining TCP/IP parameters from a central server.	RFC 2131 and RFC 2132
	FTP - File Transfer Protocol	A protocol that defines how files are transferred from one host to another on a network.	RFC 959
	HTTP - Hypertext Transfer Protocol	A protocol that defines how files are exchanged on the Internet	RFC 1945 and RFC 2068

* = IPv4 (current standard), ** = IPv6 (proposed new standard)

• • 145

The Network Interface Layer

The lowest layer in the TCP/IP stack is the Network Interface Layer. The primary responsibility of this layer is to define how a host device (computer, instrument, etc.) connects to the network. The Network Interface Layer acts as a host's connection (interface) to the network. There are no TCP/IP protocols associated with the Network Interface layer.

The Network Interface Layer is used to send and receive packets. At the Network Interface Layer, a header that contains addressing information is applied to each packet. A part of this header is the host **hardware address**. The hardware address must be unique to the device and does not change during the life of the device.

The hardware address, also called the Media Access Control (MAC) Address, Ethernet Address, Physical Address, or Network Interface Card (NIC) Address, is a 12-digit hexadecimal address. A typical hardware address is 00:30:D3:00:00:23, where the first six digits represent the manufacturer of the device and the last six digits represent the serial number assigned to the device.

NOTE

For a TCP/IP packet to be delivered to a device, the packet **must** contain the destination device's hardware address.

As packets are sent through the network, each host on the network looks at the packet to see if the packet is addressed to the host's hardware address. If not, the host ignores the packet.



The Internet Layer

The Internet Layer of the TCP/IP model contains the protocols responsible for addressing and routing of packets. The Internet Layer includes several protocols, including:

- Internet Protocol (IP)
- Address Resolution Protocol (ARP)
- Internet Control Message Protocol (ICMP)

For TCP/IP communications to be successful, the packet examined by the Network Interface Layer must include a hardware address. As the packet moves up to the Internet Layer, it also must include an **IP address**. The Internet Layer provides the protocols to determine the hardware address for routing the packet to its destination.

Internet Protocol (IP) is responsible for determining the source and (IP) Internet Protocol (IP) is responsible for determining the source and destination IP addresses of every packet sent on the network. Typically, IP addresses are assigned by a **System Administrator**. An IP address is assigned to each host on a network and each host must have a unique IP address.

In contrast to a hardware address that refers to a physical network interface card, an IP address is a **logical address** that is assigned to the host. The IP address can be changed, since it refers only to the host. A typical IP address is 167.155.21.45. See <u>IP Addressing</u> for a description of IP addresses.

- Address Resolution Protocol (ARP) Address Resolution Protocol (ARP) is a protocol used to **resolve** (translate) a **logical address** to a physical (hardware) address. ARP is used when a source host wants to communicate with a destination host, but has only the IP address. After the hardware address is resolved, ARP maintains that information for a short period of time.
- Internet Control Internet Control Message Protocol (ICMP) is primarily used to send error messages, perform diagnostics, and control data flow. You can also use ICMP with the Packet InterNet Groper (**Ping**) utility. To ping another host from a Command Prompt, type **ping** *ip* address>.



Example: Using ping for Echo Request

For example, you could use the Ping utility to send ICMP echo request packets to the destination host and request the destination host return these packets. If the packets are returned, you can assume the connection is good. If the packets are not returned, a connectivity problem exists. This figure shows an example return from pinging a computer at IP address 156.140.72.1.

```
C:>>ping 156.140.72.1

Pinging 156.140.72.1 with 32 bytes of data:

Reply from 156.140.72.1: bytes=32 time<10ms TTL=254

Ping statistics for 156.140.72.1:

Packets: Sent = 4, Received = 4, Lost = 0 (0% loss),

Approximate round trip times in milli-seconds:

Minimum = 0ms, Maximum = 0ms, Average = 0ms

C:>>
```

The Transport Layer

The Transport Layer of the TCP/IP model contains the protocols responsible for addressing and routing of packets. The Transport Layer determines if the sender and receiver will establish a connection before communicating and how often acknowledgements of the connection are sent. The protocols at the Transport Layer deliver data to and receive data from the Transport Layer protocols of other hosts. Hosts can be on the same local network or on a remote network. The Transport Layer has two protocols:

- Transmission Control Protocol (TCP)
- User Datagram Protocol (UDP)

Transmission Control Protocol (TCP) Transmission Control Protocol (TCP) is the protocol that connects the sending host and the receiving host to each other. TCP provides the connection as packets are moved between hosts having the conversation. Every packet has a TCP header that includes **sequence numbers**, acknowledgement numbers, address information, and other information. If packets get out of order enroute from the source to the destination, the sequence number allows the packets to be reassembled in the correct order.



As each packet is sent from the source host, an acknowledgement of receipt is sent by the destination host within a specified time. If the acknowledgement is not sent within this time, the sender re-sends the packet.

If the receiver gets the packet in a damaged condition, the packet is discarded and the receiver sends no acknowledgement. In this case, since an acknowledgement was not received in the specified time, the sender resends the packet. This figure shows a network in which the source and destination are using TCP and summarizes the steps to use TCP for communication between hosts.



User Datagram Protocol (UDP)

User Datagram Protocol (UDP) is used for connectionless, non-guaranteed communication. Unlike TDP, UDP does not set up a connection and does not use acknowledgements. Instead, UDP sends out packets without waiting for acknowledgement from the destination.



The Application Layer

The Application Layer is the layer where requests for data or services are processed. Applications at this layer are waiting for requests to process and all applications are "listening" at their respective **port**. The Application Layer has two protocols:

- File Transfer Protocol (FTP)
- Hypertext Transfer Protocol (HTTP)

NOTE

The Application Layer is not where an Internet browser, spreadsheet, etc. interact. Applications running at the Application Layer interact with the browser, spreadsheet, etc. applications.

Ports

This figure shows the TCP/IP protocol stack. In the stack, TCP and UDP each have access to 65,536 ports at the Application Layer. As a packet moves up the stack, IP directs the packet to either a TCP port or to a UDP port. Since all applications are listening at their respective ports, when the packet arrives at the appropriate port, it can be processed correctly.





Sockets For TCP/IP protocol, when a packet is delivered to a specific IP address, it is passed up to TCP or UDP and then to the appropriate host. This process forms a funnel through the TCP/IP stack, called a **socket**. A socket is uniquely defined by the IP address, the end-to-end protocol (TCP or UDP), and the port number.

This figure shows typical socket architecture. When a socket is first created, it has an associated protocol but not an IP address or port number. A socket must be **bound** to a port number before it can receive messages from a remote application.



• • 151

File Transfer Protocol (FTP) is an application used to transfer files from a host to another host and store the files on the requesting host. In an FTP session, one host (the client) requests a file and the other host (the server) transfers a copy of the file to the client. The two hosts thus establish a **client/server** relationship. File transfer can be in text or binary format.

An FTP client can use a word processor, FTP command line utility or FTP command server to request a file from the server. The command-line FTP client application requires you to know FTP commands, while an FTP command interpreter allows you to connect to an FTP server without using FTP commands. You can login to FTP using anonymous so that no password is required.

Hypertext Transfer Protocol (HTTP) Hypertext Transfer Protocol (HTTP) is an application used to transfer files from one host (the server) to another host (the client) and to display the files at the requesting host. The HTTP application runs on a Web server, listens for requests at a TCP port (usually port 80 for requests), and sends files back to the requestor.

The requesting host displays the files on a web browser, such as Internet Explorer. A client makes an HTTP request by typing a **Uniform Resource Locator (URL)** in the address line of the web browser or by clicking a **hyperlink** on a page that is displayed on the web browser.



IP Addressing

This section describes IP addressing for TCP/IP networks, including:

- IP Address Classes
- Subnets and Subnet Masks
- Local and Remote Networks
- IP Address Configuration Methods

IP Address Classes

Each host on a TCP/IP network must have a unique address. This address is called the **IP address** and consists of a **network portion** and a **host portion**. The network portion and host portion of an IP address are determined by the **subnet mask**.

Each IP address consists of four decimal numbers separated by periods. However, TCP/IP uses the 32-bit binary equivalent of the IP address. For example, the decimal value and binary equivalent of IP address 14.230.26.116 is:

14.230.26.116 = 0000 0110.1110 0110.0001 1010.0111 0100

IP addresses are divided into five classes: Class A, Class B, Class C, Class D, and Class E. This table summarizes standard Class A, Class B, and Class C IP addresses. In the table, the **Range** is the range of values for the first 8 bits of the IP address, regardless of class. In the **Binary Equivalent** column, $\mathbf{N} = \mathbf{a}$ network bit and $\mathbf{H} = \mathbf{a}$ host bit.

Class	Range	Network Bits	Binary Equivalent
А	1 - 127	8	NNNN NNNN.НННН НННН.НННН НННН.НННН НННН
В	128 - 191	16	NNNN NNNN.NNNN NNNN.HHHH HHHH.HHHH HHHH
С	192 - 223	24	NNNN NNNN.NNNN NNNN.NNNN NNNN.HHHH HHHH



NOTE

Dot-notation addresses ("nnn.nnn.nnn" where "nnn" is a byte value) such as IP addresses must be expressed with care, as most web software on the PC will interpret byte values with leading zeros as octal numbers. Thus, "255.255.020.011" is actually equivalent to decimal "255.255.16.9" rather than "255.255.20.11" because ".020" is interpreted as "16" expressed in octal and ".011" as "9". To avoid confusion, use only decimal expressions of byte values (0 to 255), with **no leading zeros**.

For example, the Agilent 33220A assumes that all IP addresses and other dot-notation addresses are expressed as decimal byte values and strips all leading zeros from these byte values. Thus, trying to enter "255.255.020.011" in the IP address field, it becomes "255.255.20.11" (a purely decimal expression). You should enter exactly the same expression, "255.255.20.11" in your PC web software to address the instrument. Do not use "255.255.020.011" as the PC will interpret that address differently due to the leading zeros.

Example: Class A, B, and C IP Addresses

This table shows three example IP addresses, where the bold part of the binary equivalent is the network portion of the address and the non-bold portion is the host portion of the address. For a Class A address, the first 8 bits are the network portion and the remaining 24 bits are the host portion of the address.

For a Class B IP address, the first 16 bits are the network portion and the remaining 16 bits are the host portion of the address. For a Class C IP address, the first 24 bits are the network portion and the remaining 8 bits are the host portion of the address.

Class	IP Address	Binary Equivalent
А	54.16.23.1	0011 0110 . 0001 0000. 0001 0111. 0000 0001
В	154.16.23.1	1001 0110. 0001 0000 . 0001 0111. 0000 0001
С	204.16.23.1	1100 1100. 0001 0000. 0001 0111. 0000 0001



Subnets and Subnet Masks

As noted, an IP address consists of a network portion and a host portion. A **subnet mask** is a number that looks like an IP address that shows IP how many bits are used for the network portion of the IP address by "masking" the network portion of the IP address. Every IP address must have a subnet mask. You can use standard or custom subnet masks.To see how subnet masks are used, we will first define a **subnet**.

- What is a Subnet? Large networks can be divided by a System administrator by creating a number of smaller networks, called **subnets**. The larger network is divided into subnets by using **routers** (also called **default gateways**). All devices on one side of the router form a separate subnet. Each subnet must have a unique set of IP addresses that are defined by the **subnet mask** for that network. Each side of the router must have a separate network address.
- Standard SubnetA subnet mask (like an IP address) is 32 bits long, with several contiguousMasksbits (all 1s) that represent the network portion of the IP address and the rest
of the bits (all 0s) that represent the host portion of the IP address.

For example, the standard subnet mask for a Class A address is **255**.0.0.0 (binary **1111 1111**.0000 0000.0000 0000.0000 0000), so the first eight bits (all 1s) are used for the network portion of the IP address and the remaining 24 bits (all 0s) are used for the host portion of the IP address. This table shows the standard subnet masks for Class A, Class B, and Class C IP addresses.

Class	Standard Subnet Mask	Standard Subnet Mask (Binary Equivalent)
А	255 .0.0.0	1111 1111. 0000 0000. 0000 0000. 0000 0000
В	255.255 .0.0	1111 1111. 1111 1111 . 0000 0000. 0000 0000
С	255.255.255 .0	1111 1111. 1111 1111. 1111 1111 . 0000 0000



Custom Subnet Masks You can create custom subnet masks for Class A, Class B, and Class C IP addresses. For example, for a Class A IP address, by using custom subnet mask **255.255.224**.0 you can create 2,046 unique networks with 4,094 unique host on each network.

A custom subnet mask adds more bits for the network portion of the IP address and uses fewer bits for the host portion of the IP address. These additional bits are called **subnet bits**.

NOTE

Consult a standard network or TCP/IP book or use a subnet calculator (available on the Web) for creating custom subnets.

Local and Remote Networks

For every packet sent from a source host, IP must determine if the destination host is on the same **local network** or is on a **remote network**. If the destination host is on the same local network as the source host, IP directly gets the hardware address of the destination host. If the destination host is on a remote network, IP sends packets to the **router**.

Example: Creating Local and Remote Networks

For example, consider the two networks (subnets) in this figure that are created using a router. Computer A, instrument B, instrument C, and one side of the router form one local network, while computer D, instrument E, instrument F, and the other side of the router form another local network.

If computer A wants to send a packet to instrument B or to instrument C, the destination is **local** (all hosts are on the same side of the router). In this case, IP can get the hardware address of the destination host (instrument B or instrument C) and can establish communication directly with the destination host.

However, if computer A wants to send a packet to instrument E or to instrument F, the destination is **remote** (hosts are separated by the router). In this case, IP gets the hardware address of the router. The packet is then sent to the router for transmission to the destination host (instrument E or instrument F).







IP Address Configuration Methods

This section introduces some IP address configuration methods, including:

- Configuration Methods Overview
- Dynamic Host Configuration Protocol (DHCP)
- Auto-IP/ZEROCONF
- Duplicate IP Address Detection

Configuration Methods Overview

IP address configuration methods can be divided into two categories: automatic IP configuration and manual IP configuration.

- Dynamic Host Configuration Protocol (DHCP) provides automatic TCP/IP setup for devices on networks that support DHCP.
- Auto-IP/ZEROCONF provides automatic TCP/IP setup for devices on networks that do not support DHCP.
- Manual IP Configuration provides manual TCP/IP setup for devices on manually configured networks.

This table shows the applicablility of these methods to various network topologies.

Network Topology		Manual IP		
	DHCP	Auto-IP	ZEROCONF	
Site LAN [with DHCP]	Yes	No	No	Yes
Site LAN [without DHCP]	No	Yes	Yes	Yes
Private LAN [with DHCP]	Yes	No	No	Yes
Private LAN [without DHCP]	No	Yes	Yes	Yes



Dynamic Host Configuration Protocol (DHCP)

Dynamic Host Configuration Protocol (DHCP) allows each host on the network to be automatically assigned a unique IP address when the device is connected to a network that supports DHCP and is turned ON.

To enable DHCP on a Site LAN, the System Administrator assigns a pool of IP addresses to a **DHCP Server** to be leased to hosts on the network. Each host on the network is called a **DHCP Client**.

Auto-IP/ZEROCONF

Auto-IP/ZEROCONF can be used to automatically assign IP addresses on networks that do not have DHCP Servers. Auto-IP is a defacto standard for IP address allocation implemented by Microsoft and Apple. ZEROCONF is an Internet standard (RFC) that is currently under development and is expected to be widely implemented in Linux, Windows, etc. Auto-IP and ZEROCONF are similar techniques for IP address allocation and can interoperate together.

Both Auto-IP and ZEROCONF allocate IP addresses from the link-local IP address range (169.254.xxx.xxx). The addresses are allocated using an ARP-based protocol to determine is the desired address is already in use. Auto-IP and ZEROCONF can both be used such that DHCP is tried first and then Auto-IP/ZEROCONF are used after DHCP times out (typically after about 45 - 60 seconds).

Check your instrument's *User's Guide* to determine if the instrument supports Auto-IP/ZEROCONF and, if supported, if Auto-IP/ZEROCONF can be disabled.

Duplicate IP Address Detection

Duplicate IP address detection ensures that an instrument does not attempt to start operation using an IP address that is already in use on that network. Duplicate IP Address detection provides basic diagnostic information to identify a problem on the network.

Duplicate IP Address detection is described by several Internet standards (RFC 2642/RFC 2131/ZEROCONF). It is implemented using ARP (Address Resolution Protocol) and uses broadcast Ethernet communication. The scope of duplicate IP address detection is the current Ethernet subnet. Most operating systems (Windows, Apple, etc.) implement duplicate IP address detection.



Duplicate IP Addresses on non-DHCP Networks Generally, duplicate IP addresses only occur on a manually configured IP address network that does not use DHCP and Auto-IP/ZEROCONF. For example, a user may try to determine an IP address for a host by PINGing IP addresses on the network until an IP address is selected that does not respond.

Based on this non-response, the user may assume the IP address is unused and assign it to their device. However, at a later date, a host that legitimately has that IP address will attempt to use it. As a result, both hosts will experience network problems that can be very difficult to track down.

A duplicate IP address can also happen when a user makes an error during manual entry of an IP address and accidentally configures an IP address already in use elsewhere.

Duplicate IP Addresses on DHCP Networks

Duplicate IP addresses on DHCP systems are unlikely but they are possible. The DHCP specification (RFC 2131) specifies how a duplicate IP address check should be done within the DHCP Discover/Offer/Request/ Acknowledgement protocol sequence. Ideally, the network interface should be disabled when a duplicate IP address is detected. Windows informs the user of the duplicate IP address, but does NOT disable the network interface.

If two or more devices on the network have the same IP address, the first device that starts using the duplicate IP address will not detect a problem because it is the only device using that IP address. However, when a second device starts and attempts to use the duplicate IP address, a duplicate IP address error is reported.

To recover from a duplicate IP address condition, eliminate the duplicate IP address condition and then restart the affected devices. For example, an instrument can be restarted (to clear up the error) by disconnecting the Ethernet port long enough to be detected by the Ethernet connection monitoring or by cycling power to the instrument.



Device Hostname Services

This section introduces some device Hostname services, including:

- Device Hostname Services Overview
- Dynamic DNS Naming
- RFC NetBIOS Naming
- Static DNS Naming

Device Hostname Services Overview

For ease of operation and communication on a TCP/IP network, you can assign a name to a host, called a **hostname**. However, since TCP/IP must use IP addresses to connect a host to another host, hostnames must be resolved (translated) to an IP address before TCP/IP can use the names.

Device hostname services can be divided into two categories: automatic hostname services and manual hostname services. This table shows the applicablility of these services to various network topologies.

	Automatic Hostname Se	Manual Hostname	
Network Topology	Dynamic DNS	RFC NetBIOS	Static DNS
Site LAN (Direct Connection)	Works with networks that have Dynamic DNS	Works in all network topologies	Works in all network topologies
Site LAN (Router Connection)	Works with networks that have Dynamic DNS	Works in all network topologies	Works in all network topologies
Private LAN (Router Connection)	Works with networks that have Dynamic DNS	Works in all network topologies	Works in all network topologies
Private LAN (Direct Connection)	Typically not used in this topology	Works in all network topologies	Works in all network topologies



Dynamic DNS Naming

Dynamic Domain Name System (Dynamic DNS) is a distributed database of **hostnames** and associated IP addresses on the Internet. All hostnames on the Internet are divided into categories, called **domains**, such as .com, .edu, .org, etc. Dynamic DNS automatically provides hostnames and domain names for devices on networks that support Dynamic DNS. This figure shows a typical process for a DHCP Client to register a hostname and a PTR (pointer) name via a DHCP Server.

If the registration process fails, the DHCP Client hostname may not be assigned to the instrument. If you want to change the hostname after the registration process initiated by the DHCP server, you will probably need to cycle power on the instrument to re initiate a DHCP request and hostname registration attempt.





RFC NetBIOS Naming

RFC NetBIOS Naming is a peer-to-peer naming protocol used by Microsoft File/Print Sharing that automatically provides hostnames for devices on networks that support RFC NetBIOS Naming. RFC NetBIOS naming uses a six-step process to resolve an IP address for a specified host name, as shown in the following figure. TCP/IP completes each step in the sequence shown before returning an error message. RFC NetBIOS naming does require a reboot to force a naming change to take effect.



◀ ◀

163

Static DNS Naming

In contrast to automatic DNS name resolution, **Static DNS naming** does not require any host functionality to support dynamic methods of hostname resolution. Static DNS Naming uses a seven-step process to resolve an IP address to a specified hostname, as shown in this figure. TCP/IP completes each step in the sequence shown before returning an error message. See your instrument's *User's Guide* to determine if your instrument supports Static DNS Naming.





Configuring Your PC for LAN Operation

This section shows steps to configure your PC, as required, for operation on a Private LAN or for operation on site LANs that do not support DHCP, including:

- Checking PC Settings
- Installing Network Interface Cards
- Installing TCP/IP on Your PC
- Setting PC IP Address
- Setting PC Hostname

NOTE

If your PC is part of an existing private LAN or site LAN, you probably do not need to do the steps in this section. If local policies require it, contact your System Administrator before changing PC settings.

Checking PC Settings

To check the current settings for your PC, type **ipconfig /all** from the Command Prompt. This figure shows a typical display for a Windows 2000 PC. The display for other operating systems may be different. If the settings are correct for your LAN operation, do not do any of the steps in this section. If the settings need to be changed, go to Installing Network Interface Cards.





Installing Network Interface Cards

Network Interface Cards (NIC) provide the hardware interface between your PC and network devices such as routers, hubs, or switches. As required, install NIC(s) in your PC according to the NIC manufacturer's instructions.

Installing TCP/IP on Your PC

To configure a network interface card, the TCP/IP protocol must be installed and configured. See the following examples for a Windows 2000 PC or for a Windows XP PC. Modify the steps as required for your operating system.

Example: Installing TCP/IP (Windows 2000)

- 1 Click Start | Settings | Network and Dial-up Connections to display the Network and Dial-up Connections dialog box.
- 2 Right-click Local Area Connection and then click Properties to display the Local Area Connection Properties dialog box.
- 3 The General tab should display Internet Protocol (TCP/IP). If not, click Install..., then select Protocol and click Add. Then, select TCP/IP Protocol and click Install.

ocal Area Connection	Properties	?	×
General			
Connect using:			
💷 Realtek RTL81	39(A)-based PCI Fast E	thernet Adapter	
		Configure	
Components checked	are used by this conne	ection:	
🗹 🔜 Client for Micro	osoft Networks		1
File and Printe File and Printe File and Printe	r Sharing for Microsoft col (TCP/IP)	Networks	
Install	Uninstall	Properties]
Description			
Transmission Contro wide area network across diverse inter	ol Protocol/Internet Pro protocol that provides o connected networks.	tocol. The default communication	
Show icon in task	bar when connected		1
			_



Example: Installing TCP/IP (Windows XP)

- 1 Click Start | Network| Internet Connections. From the "or pick a control panel icon", select Network Connections.
- 2 Right-click Local Area Connection and then click Properties to display the Local Area Connection Properties dialog box.
- 3 The General tab should display Internet Protocol (TCP/IP). If not, click Install, then select Protocol and click Add. Then, select TCP/IP Protocol and click Install.

General	Authentication Advanced
Connec	t using:
119 3	Com 3C920 Integrated Fast Ethernet Controller (3C905C-
	Configure
This c <u>o</u>	nnection uses the following items:
	Client for Microsoft Networks
	File and Printer Sharing for Microsoft Networks
	Juos Packet Scheduler
	Intellict Process (Per 21.)
	nstall Uninstall Properties
- Desc	iption
Tran	smission Control Protocol/Internet Protocol. The default
wide	area network protocol that provides communication ss diverse interconnected networks
📃 Sho	w icon in notification area when connected



Setting PC IP Address

Depending on the LAN capabilities, you can select automatic or manual methods to set the IP address and the DNS Server address on your PC. These examples show how to set an IP Address for a Windows 2000 or Windows XP PC. Modify the steps as required for your operating system.

NOTE

You should **not** change the IP configuration of your PC unless you are sure this is necessary. In almost all cases, the PC should be automatically configured.

Example: Setting PC IP Address (Windows 2000)

- 1 Click Start | Settings | Network and Dial-up Connections to display the Network and Dial-up Connections dialog box.
- 2 Right-click Local Area Connection and then click Properties to display the Local Area Connection Properties dialog box.
- From the General tab, select TCP/IP Protocol and then click Properties to display the Internet Protocol (TCP/IP) Properties dialog box. Select the method to set the IP address and/or Preferred and Alternate DNS server addresses.
- 4 If you select manual assignments, enter the desired values. When you have finished, click the **OK** box to enable the assignments. This figure shows an **Internet Protocol (TCP/IP) Properties** dialog box to manually assign IP and DNS Server configuration on a Windows 2000 PC.



168

Example: Setting PC IP Address (Windows XP)

- 1 Click Start | Network| Internet Connections.
- 2 From the "or pick a control panel icon", select **Network Connections**.
- 3 Right-click Local Area Connection and then click Properties to display the Local Area Connection Properties dialog box
- 4 From the General tab, select TCP/IP Protocol and then click Properties.
- 5 Select **Use the Following IP Address** and then enter the IP address and subnet mask for the network. As required, enter the default gateway and DNS settings. Follow the on-screen instructions and reboot your PC as required.
- 6 This figure shows an Internet Protocol (TCP/IP) Properties dialog box to automatically obtain IP addresses and to manually assign DNS Server addresses on a Windows XP PC.

For this PC, the IP address is automatically assigned by the network.	Internet Protocol (TCP/IP) Properties
For this PC, the preferred and alternate DNS server IP addresses are manually set.	Detain DNS server address automatically Usg the following DNS server addresses: Preferred DNS server: 156 . 140 . 113 . 8 Alternate DNS server: 156 . 140 . 2 . 12 Advanced DK Cancel

I

Setting PC Hostname

As desired, you can set an approved Hostname and a Domain Name (such as *dept.company.com*) for your PC. Examples for Windows 2000 and Windows XP follow. Modify the steps as required for your operating system.

Example: Setting PC Hostname (Windows 2000)

- 1 Right-click the **My Computer** icon and then click **Properties** to display the **System Properties** dialog box. Select the **Network Identification** tab and then click **Properties** to display the **Identification Changes** dialog box.
- 2 As required, set/change the Hostname and/or Domain Name. When you have finished, click the **OK** box to enable the assignments. This figure shows an example **Identification Changes** dialog box on a Windows 2000 PC. As desired, you can change the Hostname and/or Domain Name.

dentification Changes	?)
You can change the name and the members computer. Changes may affect access to ne	ship of this twork resources.
Computer name:	
a-04401a-120	
Full computer name: a-34401a-123.LVLD.COM	
	More
Member of	
Oomain:	
Workgroup:	
J	
OK	Cancel



Example: Setting PC Hostname (Windows XP)

- 1 Right-click the **My Computer** icon and then click **Properties** to display the **System Properties** dialog box.
- 2 Select the Network Identification tab and then click Properties to display the Identification Changes dialog box.
- 3 As required, set/change the Hostname and/or Domain Name. When you have finished, click the **OK** box to enable the assignments.
- 4 This figure shows an example **Identification Changes** dialog box on a Windows XP PC. As desired, you can change the Hostname and/or Domain Name.

System R	estore	Autom	atic Updates	Remote
General	Com	puter Name	Hardware	Advanced
Wir on J	ndows uses the network	the following in	formation to ide	ntify your computer
Computer <u>d</u> es	cription:	mypc1		
		For example: " Computer".	Kitchen Compu	ter" or "Mary's
ull computer	name:	mypc1.myname	AGILENT.COM	4
)omain:		LPC.AGILENT.	COM	
υ. Γο rename thi	s computer	or join a domair	n, click Change.	Change
			Ca	ncel Apply





Guide Information

- Guide Contents
- Related Documentation
- Accessing an Electronic Copy of This Guide
- General Information for This Guide
- Contacting Agilent



Guide Contents

This Agilent Technologies USB/LAN/GPIB Interfaces Connectivity Guide shows how to connect instruments to USB, LAN, and GPIB interfaces and how to configure and troubleshoot these interfaces on PCs with Windows 98, Windows Me, Windows NT 4.0, Windows 2000, or Windows XP operating systems. A summary of the guide contents follows.

NOTE

This guide does not describe LAN networks that include a gateway, such as the Agilent E5810A LAN/GPIB Gateway for Windows. See the applicable gateway documentation for information on gateway systems.

Also, this guide does not provide a detailed description of LAN, USB, or GPIB interfaces or TCP/IP networks. Consult standard reference texts for this information.

Section	Description			
CONNECTIVITY GUIDELINES				
Getting Started	Shows steps to get started and how to contact Agilent			
Connecting Instruments to LANs	Gives guidelines to connect instruments to LANs			
Connecting Instruments to USB	Gives guidelines to connect instruments to USB			
Connecting Instruments to GPIB	Gives guidelines to connect instruments to GPIB			
Programming Your Instruments	Gives guidelines to program instruments via LAN/USB/GPIB			
REFERENCE INFORMATION				
Installing I/O Software	Shows how to install I/O software on your PC			
Troubleshooting Guidelines	Troubleshooting guidelines for LAN, USB, and GPIB Interfaces			
TCP/IP Network Basics	Provides an overview of TCP/IP networks			
Guide Information	Lists general information for this guide			
Glossary	Defines some of the terms used in this guide			
Index	Index of selected items in this guide			



Related Documentation

Suggested related documentation you can use for interface connectivity operation follows. After the Agilent IO Libraries have been installed on your PC, .pdf files of the *Agilent IO Libraries*, *VISA User's Guide*, and *SICL User's Guide* are available. Click the blue IO icon on the Windows taskbar and then click **View Documentation**.

Product	Related Documentation
Agilent IO Libraries	For additional information on the Agilent IO Libraries, see the Agilent IO Libraries Installation and Configuration Guide for Windows.
VISA	To develop and use VISA applications, see the Agilent VISA User's Guide.
SICL	For SICL I/O applications, see the Agilent SICL User's Guide for Windows.
Programmer's Toolkit	To develop programs in Visual Studio.Net, see the <i>Agilent T&M Programmer's Toolkit for Visual Studio.NET</i>

Accessing an Electronic Copy of This Guide

- Access from the IO After the Agilent IO Libraries are installed, a blue IO icon appears on the right hand side of the Windows toolbar. To access an electronic version of this connectivity guide, click the IO icon, then click **View Documentation** and then click **Connectivity Guide**.
- Access from the CD You can access an electronic copy of this guide from the *Documentation* and Software (or equivalent title) CD that shipped with your instrument. To access an electronic version of this guide, insert the CD into your CD-ROM drive. Then, from Windows Explorer, navigate to <*drive>*/Manuals (where <*drive>* is your CD-ROM drive) and click connectivity.pdf.
- Access from the IO Libraries CD You can also access an electronic copy of this guide from the Agilent IO Libraries for Instrument Control CD. Insert the Agilent IO Libraries for Instrument Control CD into your CD-ROM drive. When the InstallShield Wizard appears, quickly click **Cancel**. Then, navigate to <**drive**/**Manuals** (where <**drive**> is your CD-ROM drive) and click **connectivity.pdf**.



General Information for This Guide

Notice

The information contained in this document is subject to change without notice.

Agilent Technologies shall not be liable for any errors contained in this document. *Agilent Technologies makes no warranties of any kind with regard to this document, whether express or implied. Agilent Technologies specifically disclaims the implied warranties of merchantability and fitness for a particular purpose.* Agilent Technologies shall not be liable for any direct, indirect, special, incidental, or consequential damages, whether based on contract, tort, or any other legal theory, in connection with the furnishing of this document or the use of the information in this document.

- U.S. Government Restricted Rights If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.
- Trademark Information Adobe[®] and Acrobat[®] are trademarks of Adobe Systems Incorporated. Microsoft[®], Visual Basic[®], Windows[®] 98, Windows[®] Me, Windows NT[®], Windows[®] 2000 and Windows[®] XP are U.S. registered trademarks of Microsoft Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

CopyrightThis Agilent Technologies USB/LAN/GPIB Connectivity Guide is © AgilentInformationTechnologies, Inc. 2003.



Contacting Agilent

You can reach Agilent Technologies at these telephone numbers:

Americas Call Center:	1-800-452-4844
Canada Call Center:	1-877-894-4414
European Call Center:	+31-20-547-9900
Japan Call Center:	+81-426-56-7832

For other countries, contact your country's Agilent support organization. A list of contact information for other countries is available on the Agilent Internet site: www.agilent.com/find/assist A list of other Agilent Websites follows.

URL	Description
www.agilent.com/find/assist	Agilent Technologies "Contact Us" page
www.agilent.com/find/iolib	Update the Agilent IO Libraries software
www.agilent.com/find/iolib_support	The latest customer support information
www.agilent.com/find/ADN	Connectivity resources all in one place
www.agilent.com/find/techsupport	Technical support information, including manuals, application notes, FAQs, and software and firmware downloads
www.agilent.com/find/connectivity	For connection, communication and control of test instruments from your computer, you can find out the latest in connectivity.





Glossary



Α	address A string uniquely identifying a particular interface or a device on that interface to identify the interface or device.
В	Bridge In telecommunication networks, a bridge is a product that connects a local area network (LAN) to another local area network that uses the same protocol (for example, Ethernet or token ring). You can envision a bridge as being a device that decides whether a message from you to someone else is going to the local area network in your building or to someone on the local area network in the building across the street. A bridge examines each message on a LAN, "passing" those known to be within the same LAN, and forwarding those known to be on the other interconnected LAN (or LANs).
С	client Part of the client/server model used in distributed computing. A client is a computer system that requests services from a server computer system, such as I/O application requests, networking requests, etc.
	controller A computer used to communicate with a remote device such as an instrument. In the communications between the controller and the device, the controller is in charge of and controls the flow of communication that is, it does the addressing and/or other bus management).
D	device A unit that receives commands from a controller. Typically a device is an instrument but could also be a computer acting in a non-controller role, or another peripheral such as a printer or plotter.
	DHCP Short for Dynamic Host Configuration Protocol, a protocol for assigning dynamic IP addresses to devices on a network. With dynamic addressing, a device can have a different IP address every time it connects to the network.
	In some systems, the device's IP address can even change while it is still connected. DHCP also supports a mix of static and dynamic IP addresses.

┥ ┥ 180 🕨
Dynamic addressing simplifies network administration because the software keeps track of IP addresses rather than requiring an administrator to manage the task. This means that a new computer can be added to a network without manually assigning it a unique IP address. DHCP client support is built into Windows 98 and NT workstations.

DNS

Е

G

Short for Domain Name System (or Service), an Internet service that translates domain names into IP addresses. Because domain names are alphabetic, they are easier to remember. The Internet, however, is really based on IP addresses. Every time you use a domain name, therefore, a DNS service must translate the name into the corresponding IP address.

For example, the domain name www.example.com might translate to 198.105.232.4. The DNS system is a distributed system. If one DNS server does not know how to translate a particular domain name, it asks another one, and so on, until the correct IP address is returned.

Ethernet (MAC) Address

The Media Access Control Address, also known as the link-level address, the Ethernet (station) Address, the LANIC ID and Hardware Address. This is a unique 48-bit address assigned by the manufacturer for each Ethernet device. It is usually displayed as 12 hexadecimal characters, sometimes with colon or dash separators between every two characters, such as "00:03:d3:00:00:17" or "00-03-d3-00-00-17".

gateway

Hardware that permits a network connection between the LAN that your computer understands and the instrument specific interface that your device understands.

Gateway IP Address

This parameter is the IP Address of the default subnet gateway that allows an instrument to communicate with systems that are not on the local subnet. Thus, this is the default subnet gateway where packets are sent that are destined for a device not on the local subnet, as determined by the subnet mask setting. Only one default subnet gateway can be configured. A value of 0.0.0.0 indicates that no subnetting is to be done.



Hub

н

Т

A common connection point for devices in a network. Hubs are commonly used to connect segments of a LAN. A hub contains multiple ports. When a packet arrives at one port, it is copied to the other ports so that all segments of the LAN can see all packets.

- **Passive hubs** serve simply as a conduit for the data, enabling it to go from one device (or segment) to another.

- Intelligent hubs include additional features that enables an administrator to monitor the traffic passing through the hub and to configure each port in the hub. Intelligent hubs are also called manageable hubs.

- **Switching hubs** actually read the destination address of each packet and then forward the packet to the correct port.

instrument

A device that accepts commands and performs a test or measurement function.

interface

A connection and communication media between devices and controllers, including mechanical, electrical, and protocol connections.

IP Address

An Internet Protocol (IP) address is an identifier for a computer or device on a TCP/IP network. Networks using the TCP/IP protocol route messages based on the IP address of the destination. The format of an IP address is a 32-bit numeric address written as four numbers separated by periods. Each number can be zero to 255. For example, 1.160.10.240 could be an IP address.

Within an isolated network, you can assign IP addresses at random as long as each one is unique. However, connecting a private network to the Internet requires using registered IP addresses (called Internet addresses) to avoid duplicates. The four numbers in an IP address are used in different ways to identify a particular network and a host on that network. The InterNIC Registration Service assigns Internet addresses from the following three classes.

- Class A supports 16 million hosts on each of 127 networks
- Class B supports 65,000 hosts on each of 16,000 networks
- Class C supports 254 hosts on each of 2 million networks

182

LAN

L

Local Area Network. A computer network that spans a relatively small area. Most LANs are confined to a single building or group of buildings. However, one LAN can be connected to other LANs over any distance via telephone lines and radio waves. A system of LANs connected in this way is called a wide-area network (WAN).

Most LANs connect workstations and personal computers. Each node (individual computer) in a LAN has its own CPU with which it executes programs, but it also is able to access data and devices anywhere on the LAN. This means that many users can share expensive devices, such as laser printers, as well as data. Users can also use the LAN to communicate with each other, by sending e-mail or engaging in chat sessions.

There are many different types of LANs Ethernets being the most common for PCs. Most Apple Macintosh networks are based on Apple's AppleTalk network system, which is built into Macintosh computers. The following characteristics differentiate one LAN from another:

topology : The geometric arrangement of devices on the network. For example, devices can be arranged in a ring or in a straight line.

protocols: The rules and encoding specifications for sending data. The protocols also determine whether the network uses a peer-to-peer or client/server architecture.

media: Devices can be connected by twisted-pair wire, coaxial cables, or fiber optic cables. Some networks do without connecting media altogether, communicating instead via radio waves.

LANs are capable of transmitting data at very fast rates, much faster than data can be transmitted over a telephone line. However, the distances are limited and there is also a limit on the number of computers that can be attached to a single LAN.

lock

A state that prohibits other users from accessing a resource, such as a device or interface.

logical unit

A logical unit is a number associated with an interface. In Agilent SICL, a logical unit uniquely identifies an interface. Each interface on the controller must have a unique logical unit.



network

A group of two or more computer systems linked together. There are many types of computer networks, including:

local-area networks (LANs) : The computers are geographically close together (that is, in the same building).

wide-area networks (WANs) : The computers are farther apart and are connected by telephone lines or radio waves.

campus-area networks (CANs): The computers are within a limited geographic area, such as a campus or military base.

metropolitan-area networks MANs): A data network designed for a town or city.

home-area networks (HANs): A network contained within a user's home that connects a person's digital devices.

In addition to these types, the following characteristics are also used to categorize different types of networks:

topology : The geometric arrangement of a computer system. Common topologies include a bus, star, and ring. See the Network topology diagrams in the Quick Reference section of Webopedia.

protocol : The protocol defines a common set of rules and signals that computers on the network use to communicate. One of the most popular protocols for LANs is called Ethernet. Another popular LAN protocol for PCs is the IBM token-ring network.

architecture : Networks can be broadly classified as using either a peerto-peer or client/server architecture.

Computers on a network are sometimes called nodes. Computers and devices that allocate resources for a network are called servers

network protocols

The protocol defines a common set of rules and signals that computers on the network use to communicate. One of the most popular protocols for LANs is called *Ethernet*. A description of the ARP, DHCP, DNS, FTP, HTTP, ICMP, and RPC protocols follows.

ARP: Short for Address Resolution Protocol, a TCP/IP protocol used to convert an IP address into a physical address (called a DLC address), such as an Ethernet address.

DHCP: See DHCP



DNS: See DNS

FTP: Abbreviation of File Transfer Protocol, the protocol used on the Internet for sending files.

HTTP: Short for HyperText Transfer Protocol, the underlying protocol used by the World Wide Web. HTTP defines how messages are formatted and transmitted, and what actions Web servers and browsers should take in response to various commands.

ICMP: Short for Internet Control Message Protocol, an extension to the Internet Protocol (IP) defined by RFC 792. ICMP supports packets containing error, control, and informational messages. The ping command, for example, uses ICMP to test an Internet connection.

RPC: Abbreviation of remote procedure call, a type of protocol that allows a program on one computer to execute a program on a server computer. Using RPC, a system developer need not develop specific procedures for the server. The client program sends a message to the server with appropriate arguments and the server returns a message containing the results of the program executed.

proxy server

A server that sits between a client application, such as a web browser, and a real server. It intercepts all requests to the real server to see if it can fulfill the requests itself. If not, it forwards the request to the real server. Proxy servers have two main purposes:

Improve Performance: Proxy servers can dramatically improve performance for groups of users, since a proxy server saves the results of all requests for a certain amount of time. Consider the case where both user X and user Y access the Web through a proxy server. First user X requests a certain Web page, which we will call Page 1. Sometime later, user Y requests the same page.

Instead of forwarding the request to the Web server where Page 1 resides, which can be a time-consuming operation, the proxy server returns the Page 1 it already fetched for user X. Since the proxy server is often on the same network as the user, this is a much faster operation. Real proxy servers support hundreds or thousands of users.

Filter Requests: Proxy servers can also be used to filter requests. For example, a company might use a proxy server to prevent its employees from accessing a specific set of Web sites.





Router

R

S

A device that connects any number of LANs. Routers use headers and a forwarding table to determine where packets go. They use Internet Control Message Protocol (ICMP) to communicate with each other and configure the best route between any two hosts. Very little filtering of data is done through routers. Routers do not care about the type of data they handle. Routers often have DHCP Server capability.

server

Part of the client/server model used in distributed computing. The server is a computer system designated to act as a main servicer of requests from other client computer systems, such as I/O application requests, networking requests, and so forth.

SICL

The Agilent Standard Instrument Control Library, which is software used for I/O application programming.

SRQ

Service Request. An asynchronous request (an interrupt) from a remote device indicating that the device requires servicing.

Subnet

A portion of a network that shares a common address component. On TCP/IP networks, subnets are defined as all devices whose IP addresses have the same prefix. For example, all devices with IP addresses that start with 100.100.100. would be part of the same subnet. Dividing a network into subnets is useful for both security and performance reasons. IP networks are divided using a subnet mask.

Subnet Mask

A mask used to determine to what subnet an IP address belongs. An IP address has two components: the network address and the host address. For example, consider the IP address 150.215.017.009. Assuming this is part of a Class B network, the first two numbers (150.215) represent the Class B network address and the second two numbers (017.009) identify a particular host on this network.

Subnetting enables a System Administrator to further divide the host part of the address into two or more subnets. In this case, a part of the host address is reserved to identify the particular subnet.



For example, the full address for 150.215.017.009 is 10010110.11010111.00010001.00001001. The Class B network part is 10010110.11010111 and the host address is 00010001.00001001. If this network is divided into 14 subnets, the first four bits of the host address (0001) are reserved for identifying the subnet.

The subnet mask is the network address plus the bits reserved for identifying the subnetwork. (By convention, the bits for the network address are all set to 1, though it would also work if the bits were set exactly as in the network address.) In this case, the subnet mask is 1111111.1111111.11110000.00000000.

It is called a mask because it can be used to identify the subnet to which an IP address belongs by performing a bitwise AND operation on the mask and the IP address. The result is the subnetwork address:

Subnet Mask 255.255.240.000 1111111.11111111.1110000.00000000 IP Address 150.215.017.009 10010110.11010111.00010001.00001001 Subnet Address 150.215.016.000 10010110.11010111.00010000.00000000

Switch

A device that filters and forwards packets between LAN segments. Switches operate at the data link layer (layer 2) and sometimes the network layer (layer 3) of the OSI Reference Model and therefore support any packet protocol. LANs that use switches to join segments are called switched LANs or, in the case of Ethernet networks, switched Ethernet LANs. A hub connects all the devices on its "ports" together.

A switch is a bit smarter, as it understands when two devices (out of four, five, eight, sixteen, or even more) want to talk to each other and gives them a switched connection.

symbolic name

A name corresponding to a single interface. This name uniquely identifies the interface on a controller or gateway. If there is more than one interface on the controller or gateway, each interface must have a unique symbolic name.



Universal Plug and Play

Universal Plug and Play (UPnP) is an open industry standard that uses Internet and Web protocols to enable devices such as PCs, peripherals, intelligent appliances, and wireless devices to be plugged into a network and automatically know about each other.

UPnP is an architecture for pervasive peer-to-peer network connectivity of PCs and intelligent devices or appliances. UPnP builds on Internet standards and technologies, such as TCP/IP, HTTP, and XML, to enable LAN devices to automatically connect with one another and work together.

With UPnP, when a user plugs a device into the network, the device will configure itself, acquire a TCP/IP address, and use a discovery protocol based on the Internet's Hypertext Transfer Protocol (HTTP) to announce its presence on the network to other devices. UPnP devices use Extensible Markup Language (XML) to establish a common language, or "protocol negotiation", to talk to each other and determine capabilities.

VEE

The Agilent Visual Engineering Environment, which is software used for I/O application programming.

VISA

The Agilent Virtual Instrument Software Architecture library, which is software used for I/O application programming.



V

ABCDEFGHIJKLMNOPQRSTUVWXYZ

Index

Α

address, 180 Address Resolution Protocol (ARP), 147 Agilent IntuiLink, 103 Agilent IO Libraries Control icon, 34, 69 Agilent websites, 177 Auto-IP/ZEROCONF, 159

В

bridge, 180

С

client, 180 Configuring the Instrument page, 32 connecting instruments to LANs, 11 controller, 180 copyright information, 176 crossover cable, 141 crossover point, 134 Current Configuration page, 31

D

default gateways, 155 device, 180 DHCP, 159, 180 DNS, 181 domain, 162 dot-notation IP addresses, 154 dynamic host configuration protocol, 159

Ε

Ethernet (MAC) Address, 181 examples Addressing GPIB Instruments, 90 Addressing LAN Instruments, 84 Assigning a USB Alias Name, 49, 87

E (cont'd)

examples (continued) Cache and Page Refresh, 116 Class A, B, and C IP Addresses, 154 Connecting a GPIB Instrument, 66 Connecting Multiple GPIB Instrs, 67 Creating Local/Remote Networks, 156 Direct Connection (Private LAN), 26 Direct Connection (Site LAN), 24 Enabling Javascript, 117 Installing Software from a CD, 106 Installing TCP/IP (Win 2000), 166 Installing TCP/IP (Win XP), 167 Setting PC Hostname (Win 2000), 170 Setting PC Hostname (Win XP), 171 Setting PC IP Address (Win 2000), 168 Setting PC IP Address (Win XP), 169 Setting Proxy Notification, 115 Sine Wave (VB and VISA-COM), 91 Sine Wave (C++ and IVI-COM), 92 Sine Wave (C++ and VISA), 97 Sine Wave (C++ and VXIpnp), 95 Switch Connections (Private LAN), 27 Switch Connections (Site LAN), 25 Using Network Protocols, 82 Using ping Command, 112 Using ping for Echo Request, 148 Using route add Command, 113 Using VISA Assistant, 39, 53, 74 Using Web Browser, 109

F

File Transfer Protocol (FTP), 152 firewall, 141 Found New Hardware, 49, 137 Friendly Name, 181



A B C D E F G H I J K L M N O P Q R S T U V W X Y Z

G

qateway, 181 gateway IP address, 181 getting started, 7, 9 glossary, 179 **GPIB** addressing instruments, 88 Agilent IO Libraries checks, 135 Bus Address, 71 cable connection guidelines, 65 check for driver files, 130 communicate with instrs, 73 configuration files, installing, 69 configure GPIB cards, 68 connect instruments, 64 crossover point, 134 definition, 56 Found New Hardware, 137 hardware checks, 126 install GPIB cards, 60 install software, 59 interrupt mode, 133 Logical Unit, 71 polling mode, 133 program your instruments, 75 quick start, 56 SICL interface name, 71 software checks, 130 supported OS, 68 System Controller, 71 troubleshooting overview, 124 VISA Interface Name, 71

guide

contents, 174 copyright information, 176 electronic copies, 175 related documentation, 175 restricted rights, 176 trademark information, 176 using, 8 warranty, 176

н

hardware address, 146 host, 141 host portion, 153 hub, 141, 182 hyperlink, 152 Hypertext Transfer Protocol, 152

L

ihint, 133 installing software hardware requirements, 102 installation steps, 103 IVI-COM drivers, 103 pre-installation checks, 102 select Agilent IO Libraries, 104 select application software, 104 software requirements, 102 VXIppp drivers, 103 instrument Welcome page, 30 instrument, definition, 182 interface, 182 Internet Ctrl Message Protocol, 147 IO Config utility, 34, 68 iopen, 84, 90 IP addresses, 141, 154, 182 IVI-COM drivers, 103

L

LAN

Auto-IP/ZEROCONF, 159 cache and page refresh, 115 checking PC Settings, 165 checking proxy settings, 114 client-server model, 142 configuring PC for LAN operation, 165 connecting instruments, 11 **DHCP**, 159 **Display Devices Present**, 37



ABCDEFGHIJKLMNOPQRSTUVWXYZ

L (cont'd)

LAN (continued) duplicate IP address detection, 159 enabling Javascript, 117 hardware architecture, 142 installing TCP/IP, 166 interface overview, 140 IP address classes, 153 IP addressing, 153 local and remote networks, 156 Network Interface Cards, install, 166 programming instruments, 40 quick start, 12 request for comments (RFC), 145 **RFC NetBIOS Naming**, 163 setting PC hostname, 170 setting PC IP address, 168 SICL-LAN protocol, 81 Static DNS Naming, 164 subnets and subnet masks, 155 TCP/IP Application Layer, 150 TCP/IP LAN Client, configuring, 34 TCP/IP Network Interface Layer, 146 TCP/IP network model, 144 TCP/IP protocols, 81, 144 TCP/IP Transport Layer, 148 troubleshooting overview, 108 typical network topologies, 140 VXI-11 protocol, 81 lock, 183 logical address, 147 logical unit, 183

Μ

Media Access Control (MAC) address, 146

Ν

navigation arrows, 8 network classes, 182 Network Interface cards, installing, 166 network portion, 153 network protocols, 184

0

operating systems, supported, 8

Ρ

packet switching, 143 ports, 150 private LAN connecting devices, 26 definition, 12, 17 network info card, 21 programming instruments, 40 protocol stack, 144 protocol suite, 144 protocols, 144 proxy server, 185

R

request for comments (RFC), 145 restricted rights, 176 RFC NetBIOS Naming, 163 routers, 141, 186 rsrcName, 86

S

servers, 141, 186 SICL, 186 SICL-LAN protocol, 81 site LAN connecting instruments, 24 definition, 12, 16 information, gathering, 19



S (cont'd)

site LAN (continued) network info card, 20 sockets, 151 SRQ, 186 Static DNS Naming, 164 subnet, 186 subnet mask, 153, 186 switch, 141, 187 symbolic name, 187

т

TCP/IP LAN Client, configuring, 34 network basics, 139 network model, 144 protocols, 81, 144 trademark information, 176 Transmission Control Protocol. 148 troubleshooting Agilent IO Libraries checks, 135 cache and page refresh, 115 check GPIB driver files, 130 check USB driver files. 122 checking proxy settings, 114 enabling Javascript, 117 GPIB crossover point, 134 GPIB hardware checks, 126 GPIB interrupt mode, 133 **GPIB** overview, 124 GPIB polling mode, 133 GPIB software checks, 130 LAN overview, 108 USB hardware checks, 120 USB overview, 118 USB software checks, 121

U

Uniform Resource Locator, 152 Universal Plug and Play, 188 USB adressing instruments, 86 assign Alias name, 48 check instrument ID, 50 check USB driver files, 122 communicate with instrs, 52 connecting instruments, 46 definition, 42 hardware checks, 120 program your instruments, 54 quick start, 42 software checks, 121 supported OS, 42 troubleshooting overview, 118 USB Devices dialog box, 86 User Datagram Protocol (UDP), 149

Using the ping Command, 112 Using the route add Command, 113 using this guide, 8

V

VEE. 188 viOpen, 84, 90 **VISA**, 188 VISA Assistant, 39, 53 VISA rsrcName, 86 viSetAttribute, 133 VXI-11 protocol, 81 VXIppp drivers, 103

W

warranty, 176 Web Browser, 30 instrument Welcome page, 30 set TCP/IP parameters, 32 TCP/IP parameters, setting, 30 websites, Agilent, 177

